

## GL-Screw Hose Connector GL 스크류 호스 연결관

• 캡과 Connector를 Set로 공급

### 기본형 PP Connector + PBT Cap, 140°C 내열 ①

Cat. No.	Model	규격	형태	Connector od	EA Price ₩
C20-83-009	1350000S	GL 14	일자형	Ø8 mm	5,500
C20-83-018	1360000S		곡형		5,800
C20-83-027	1350001S	GL 18	일자형	Ø11 mm	13,300
C20-83-036	1360001S		곡형		16,700

### 고온용 PFA Connector + PPS Cap, 250°C 내열 ②



Cat. No.	Model	규격	형태	Connector od	EA Price ₩
C20-83-109	D581-02	GL 14	일자형	Ø8.7 mm	12,900
C20-83-118	D582-02		곡형		16,400
C20-83-127	D581-04	GL 18	일자형	Ø10.5 mm	26,300
C20-83-136	D582-04		곡형		26,300

### Quick Disconnection Type

Cat. No.	Model	규격	Connector od	EA Price ₩
C20-83-209	KA45	GL 14	Ø4.5 mm ③	9,400
C20-83-218	KA90		Ø9 mm ④	9,400

### Quick Disconnection Basic Plain Type


Cat. No.	Model	규격	Connector od	EA Price ₩
C20-83-223	KA45-TA	GL 14	Ø4.5 mm	17,600
C20-83-227	KA90-TA	Adapter	Ø9 mm ⑤	17,600



## PTFE Rigid GL Connector 테프론 GL 커넥터

Cat. No.	Model	Thread	Total Length	EA Price ₩
C20-460-751	H903-14	GL 14	49 mm	260,200
C20-460-755	H903-18	GL 18		274,100
C20-460-759	H903-25	GL 25	51 mm	440,400

- 재질 : PTFE, PPS, -20 ~ 200°C 내열
- 테프론 재질로 최상의 내화학성을 가짐
- 시료는 단지 PTFE 재질과만 접촉됨
- BOLA



## GL Screw Reducing Adapter GL 스크류 확대 / 축소형 아답터

Cat. No.	Model	From Thread GL①	To Thread GL②	Length Size	EA Price ₩
A01-460-101	H904-02	GL 18	GL 14	20.0mm	70,500
A01-460-105	H904-03	GL 25	GL 14	27.5mm	110,600
A01-460-109	H904-04	GL 25	GL 18	28.3mm	84,100
A01-460-113	H904-05	GL 32	GL 18	32.5mm	118,900
A01-460-117	H904-06	GL 32	GL 25	28.0mm	106,500
A01-460-121	H904-07	GL 45	GL 25	38.0mm	164,100
A01-460-125	H904-08	GL 45	GL 32	33.0mm	165,700

- 재질 : PTFE, -20 °C to 200°C 내열
- 다른 규격의 GL 스크류를 연결하기 위한 어답터
- 테프론 재질로 최상의 내화학성을 가짐
- PTFE/ FKM Gasket 포함
- BOLA

