




Luer Adapter, 316 Stainless Steel

스테인레스 루어 어댑터

- FLL (Female Luer Lock)
- MLL (Male Luer Lock)
- Cadence Science
- MMFLL (Micro-Mate® Female Luer Lock)
- ML (Male Luer) male luer slip that fits with a female luer



Luer-to-Luer Adapter ①

Cat. No.	Model	Side 1	Side 2	EA Price ₩
C20-55-084	6454 IND	MM FLL	ML	89,600
C20-55-086	6519 IND	FLL	MLL	164,500
C20-55-088	6520 IND	MLL	MLL	164,500
C20-55-090	6521 IND	FLL	FLL	164,500

Threaded End(NPT) Adapter ②

Cat. No.	Model	Side 1	Side 2	EA Price ₩
C20-55-091	6455 IND	MM FLL	1/8-27(NPT) Thread	113,900
C20-55-094	6526 IND	MLL	1/8-27(NPT) Thread	164,500
C20-55-097	6527 IND	MLL	1/4-18(NPT) Thread	164,500
C20-55-135	6531 IND	FLL	1/8-27(NPT) Thread	165,100
C20-55-138	6532 IND	FLL	1/4-18(NPT) Thread	165,100

Luer-to-Closed End Adapter ③

Cat. No.	Model	Side 1	Side 2	EA Price ₩
S15-205-609	6200 IND	FLL	Closed End	52,300
S15-205-611	6450 IND	MM	Closed End	164,500
S15-205-613	6517 IND	MLL	Closed End	164,500
S15-205-616	6522 IND	ML	Closed End	164,500

Threaded End(UTS) Adapter ④

Cat. No.	Model	Side 1	Side 2	EA Price ₩
C20-55-112	6458 IND	MM FLL	1/4-28 Standard Thread	126,700
C20-55-118	6524 IND	MLL	1/4-28 Standard Thread	164,500
C20-55-120	6529 IND	FLL	1/4-28 Standard Thread	164,500
C20-55-098	6456 IND	FLL	10-32 Internal Female Thread	126,700
C20-55-105	6457 IND	LL	10-32 Standard Thread	91,100
C20-55-117	6523 IND	MLL	10-32 Standard Thread	164,500

Tubing Connector ⑤

Cat. No.	Model	Side 1	Hose End 사용 Tubing id	EA Price ₩
C20-55-063	6451 IND	MM FLL	1.59 ~ 2.38mm	91,100
C20-55-070	6452 IND		3.17 ~ 4.76mm	91,100
C20-55-077	6453 IND		6.35 ~ 7.93mm	91,100
C20-55-100	6510 IND	MLL	1.59 ~ 2.38mm	164,500
C20-55-103	6511 IND		3.17 ~ 4.76mm	117,300
C20-55-115	6512 IND		6.35 ~ 7.93mm	117,300
C20-55-106	6513 IND	FLL	1.59 ~ 2.38mm	117,300
C20-55-109	6514 IND		3.17 ~ 4.76mm	117,300
C20-55-114	6515 IND		6.35 ~ 7.93mm	113,900
C20-55-102	6217 IND	FLL	3.17 ~ 4.76mm	117,300
C20-55-901	6169 IND	MLL	6.35 ~ 7.93mm	24,700

3-Way Adapter ⑥

Cat. No.	Model	Side 1	Side 2	Side Arm	EA Price ₩
C20-55-194	6163 IND	FLL	MLL	FLL	133,400
C20-55-197	6134 IND	Y Manifold FLL	MLL	MLL	125,900
C20-55-203	6136 IND	Y Tubing Connector 3 x 8mm			111,400

Peek Luer Adapter 피크 루어 어댑터

- 6601 + Inert Tube-end Fitting으로 주사기와 Teflon Tubing을 연결 시킬 수 있음
- Cadence Science



Cat. No.	Model	Luer Lock	Thread	EA Price ₩
C20-55-007	3600	Male	1/4" - 28 Internal Standard Thread Female	108,700
C20-55-014	3601	Female		108,700
C20-55-021	3602	Male	1/4" - 28 Standard Thread Male	108,700
C20-55-028	3603	Female		108,700

Inert Tube-end Fitting 튜브 연결용 피팅

- 재질 : ETFE (Cone), PP (Nut)
- Pressure : 500 psi (od1/16"), 250 psi (od1/8")
- Chemglass



Cat. No.	Model	규격	Tubing od	EA Price ₩
A01-132-664	1164-N-02	1/4"-28 UNF Male Fitting ①	1/8"	5,800
A01-132-672	1164-F-02	1Inverted Cone ②	1/8"	3,400