

Luer Connector / Quick Disconnecter

루어용 튜브 연결관, PP

- 재질 : Natural Polypropylene (-13 ~ 90°C), USP Class VI, Rohs
- Female과 Male Luer Lock을 같이 분리형 연결관처럼 사용 가능
- Tommyheco



Female Luer Lock 형 연결관 ①

Cat. No.	Model	사용 Tubing id	EA Price ₩
C20-140-624	PPLF10	1/16" / 1.59 mm	2,500
C20-140-628	PPLF15	3/32" / 2.38 mm	2,500
C20-140-632	PPLF20	1/8" / 3.18 mm	2,500
C20-140-636	PPLF25	5/32" / 3.97 mm	2,500
C20-140-640	PPLF30	3/16" / 4.76 mm	2,500
C20-140-644	PPLF40	1/4" / 6.35 mm	2,500

Male Luer Lock 형 연결관 ②

Cat. No.	Model	EA Price ₩
C20-140-652	PVLM10	3,400
C20-140-656	PVLM15	3,400
C20-140-660	PVLM20	3,400
C20-140-664	PVLM25	3,400
C20-140-668	PVLM30	3,400
C20-140-672	PVLM40	3,400

Female Luer Coupler ③

Cat. No.	Model	EA Price ₩
C20-140-648	PPLFC	1,900

Female (PPLFP) & Male (PPLMP) Plug

Cat. No.	Model	EA Price ₩
C20-140-676	PPLFP ④	4,300
C20-140-678	PPLMP ⑤	2,500




1-Way Luer Stopcocks

1-Way 루어형 스탑콕크

- 재질 : 니켈 도금 황동 본체 or 316 Stainless Steel
- Tested to 45 PSI water pressure
- Cadence Science

Nickel Plated Brass ①

Cat. No.	Model	Side 1	Side 2	EA Price ₩
C20-55-126	6019 IND		Male Luer	164,500
C20-55-128	6020 IND	Female Luer Lock	Female Luer Lock	182,400
C20-55-133	6021 IND		Male Luer Lock	237,500
C20-55-139	6021HP IND	Threaded Female Luer Lock	Male Luer Lock	386,900
C20-55-142	6021TF IND	Male Luer Lock	Threaded Female Luer Lock	164,500
C20-55-145	6023 IND		Male Luer Lock	208,800
C20-55-148	6026 IND	6mm Hose End	6mm Hose End	208,800
C20-55-154	6029 IND		6mm Hose End	240,700
C20-55-157	6030 IND	Male Luer Lock	8mm Hose End	240,700



Stainless Steel ②

Cat. No.	Model	Side 1	Side 2	EA Price ₩
C20-55-136	6021SS IND	Female Luer Lock	Male Luer Lock	605,800


3-Way Luer Stopcocks

3-Way 루어형 스탑콕크

- 재질 : 니켈 도금 황동 본체 or 316 Stainless Steel
- Tested to 45 PSI water pressure
- Cadence Science

Nickel Plated Brass ①

Cat. No.	Model	Side 1	Side 2	Side Arm	EA Price ₩
C20-55-161	6001 IND		Male Luer	Right Side,	266,100
C20-55-164	6002 IND		Female Luer Lock	Female Luer Lock	266,100
C20-55-167	6010 IND	Female Luer Lock		Right Side, 3.5mm Hose End	251,900
C20-55-170	6011 IND				212,400
C20-55-176	6011HP IND	Threaded Female Luer Lock	Male Luer Lock	Right Side, Female Luer Lock	509,000
C20-55-179	6011TF IND	Male Luer Lock		Female Luer Lock	196,500
C20-55-182	6014 IND	Male Luer Lock			319,900
C20-55-185	6017 IND	Female Luer Lock		Left Side, Female Luer Lock	319,900
C20-55-191	6056 "Y" Type IND	Female Luer Lock	Male Luer Lock	Male Luer Lock	475,900



Stainless Steel ②

Cat. No.	Model	Side 1	Side 2	Side 3	EA Price ₩
C20-55-173	6011SS IND	Female Luer Lock	Male Luer Lock	Right Side, Female Luer Lock	722,000
C20-55-188	6017SS IND	Female Luer Lock	Male Luer Lock	Left Side, Female Luer Lock	708,200