


PTFE Thin Wall Joint Sleeve 테프론 슬리브 셉타

- 재질 : PTFE (Polytetrafluorethylene)
- -200 ~ +280℃
- 0.05 mm 두께로 내약품성 / 내구성이 탁월하며, Glass Joint Cone에 씌워 완벽한 Sealing으로 진공용으로도 쓰임
- Cowie



Socket, Medium Length

Cat. No.	Model	Socket NS	PK Price ₩
S12-450-523	019.110	10/19	10ea 39,500
S12-450-529	019.112	12/21	10ea 43,100
S12-450-535	019.114	14/23	10ea 16,600
S12-450-541	019.119	19/26	10ea 16,100
S12-450-547	019.124	24/29	10ea 17,800
S12-450-553	019.129	29/32	10ea 19,100
S12-450-559	019.134	34/35	10ea 28,000
S12-450-571	019.145	45/40	10ea 40,300

Socket, Full Length

Cat. No.	Model	Socket NS	PK Price ₩
S12-450-634	019.010	10/30	10ea 60,000
S12-450-639	019.012	12/30	10ea 48,300
S12-450-644	019.014	14/35	10ea 26,700
S12-450-649	019.019	19/38	10ea 26,500
S12-450-654	019.024	24/40	10ea 20,300
S12-450-659	019.029	29/42	10ea 30,700
S12-450-664	019.034	34/45	10ea 32,900
S12-450-669	019.045	45/50	10ea 50,500

PTFE Rigid Joint Sleeve 고진공용 테프론 슬리브 셉타

- 재질 : PTFE (Polytetrafluorethylene)
- -200 ~ +280℃
- 두께 : 0.5 mm (±0.05 mm)
- 고진공용, 재사용 가능
- Leakage Rate : Less than 1×10^{-4} Torr
- Cowie



Socket, Medium Length

Cat. No.	Model	Socket NS	EA Price ₩
S12-450-723	020.010	10/19	9,600
S12-450-735	020.014	14/23	9,600
S12-450-741	020.019	19/26	11,200
S12-450-747	020.024	24/29	12,400
S12-450-753	020.029	29/32	16,300
S12-450-759	020.034	34/35	20,500

Socket, Full Length

Cat. No.	Model	Socket NS	EA Price ₩
S12-450-834	020.110	10/30	9,900
S12-450-844	020.114	14/35	10,100
S12-450-849	020.119.1	19/22	11,600
S12-450-854	020.124	24/40	15,400
S12-450-859	020.129	29/42	20,800
S12-450-864	020.134	34/45	22,300
S12-450-869	020.145	45/50	37,100

PTFE Spherical Ground Joint Sleeve PTFE 테프론 구형 슬리브

Cat. No.	Model	Size	PK Price ₩
C20-460-801	H934-02	S 13	10ea 128,300
C20-460-805	H934-04	S 19	10ea 118,100
C20-460-809	H934-06	S 29	10ea 142,700
C20-460-813	H934-08	S 35	10ea 198,400
C20-460-817	H934-12	S 40	10ea 305,100
C20-460-821	H934-16	S 51	10ea 562,200
C20-460-825	H934-18	S 64	10ea 621,100

- 재질 : PTFE, -200 °C to 250°C 내화학 및 내열
- 외부에 Gripping Ring과 Sealing Ring이 있으며 고진공 씰링용으로 적합
- BOLA



PTFE Flexible Bellow Connector 테프론 플렉서블 연결관, 120°

Cat. No.	Model	Socket & Cone	길이	EA Price ₩
C20-450-523	008.01414	14/23	90 mm	79,900
C20-450-529	008.01919	19/26	100 mm	87,800
C20-450-535	008.1919	19/38	120 mm	76,100
C20-450-537	008.02424	24/29	110 mm	98,300
C20-450-541	008.2424	24/40	130 mm	104,100
C20-450-547	008.02929	29/32	115 mm	121,400
C20-450-553	008.2929	29/42	140 mm	116,600
C20-450-559	008.03434	34/35	120 mm	133,700

- PTFE (Polytetrafluorethylene)
- -200 ~ +280 °C
- 120° 까지 휘어지는 테프론 연결관
- Cowie

