

## PTFE Screw-in Tube Fitting PTFE 스크류인 튜브 피팅

- 재질 : PTFE, -200 °C to 250°C 내열
- 사용 압력 : Max 5 bar
- 테프론 재질로 최상의 내화학성을 가짐
- 시료는 단지 PTFE만 접촉 되어짐
- BOLA


Cat. No.	Model	Thread Fitting ①	Tubing	Bore Dia	Screw-in Thread ②	EA Price ₩
C20-460-011	D518-02	GL 14	od 4 mm	Ø 4 mm	NPT 1/8"	45,800
C20-460-015	D518-04				G 1/8"	56,700
C20-460-019	D518-06				NPT 1/4"	39,300
C20-460-023	D518-08				G 1/4"	60,700
C20-460-027	D518-12		od 6 mm		NPT 1/4"	37,600
C20-460-031	D518-14				G 1/4"	54,100
C20-460-035	D518-16				NPT 3/8"	77,700
C20-460-039	D518-18				G 3/8"	99,500
C20-460-043	D518-24	GL 18	od 8 mm	Ø 8 mm	NPT 1/4"	38,300
C20-460-047	D518-26				G 1/4"	58,700
C20-460-051	D518-28				NPT 3/8"	53,600
C20-460-055	D518-30				G 3/8"	95,800
C20-460-059	D518-36		od 10 mm		NPT 1/4"	38,800
C20-460-063	D518-38				G 1/4"	58,900
C20-460-067	D518-40				NPT 3/8"	54,200
C20-460-071	D518-42				G 3/8"	96,500
C20-460-079	D518-50	GL 28	od 12 mm	Ø 12 mm	G 3/8"	136,500
C20-460-083	D518-52				NPT 1/2"	135,400
C20-460-087	D518-54		od 14 mm		G 1/2"	122,300
C20-460-095	D518-62				G 1/2"	121,000
C20-460-099	D518-68		od 16 mm		NPT 1/2"	140,500
C20-460-103	D518-70				G 1/2"	128,800



## PTFE Screw-in Tube Fitting PTFE 스크류인 튜브 피팅, 분배용

Cat. No.	Model	Inlet ①	Inlet Tubing	Outlet ②	Outlet Tubing	Bore Dia	EA Price ₩
C20-460-131	D512-01	GL14 x 1개	od 4 mm	GL14 X 3개	od 4mm	Ø 6mm	457,700
C20-460-135	D512-02		od 6 mm		od 6mm		453,600
C20-460-139	D512-08		od 4 mm	GL14 X 4개	od 4mm		692,500
C20-460-143	D512-09		od 6 mm		od 6mm		770,800

- 재질 : PTFE, -200 °C to 250°C 내열
- 사용 압력 : Max 5 bar
- 테프론 재질로 최상의 내화학성을 가짐
- 시료는 단지 PTFE만 접촉 되어짐
- BOLA



## PTFE Plug for Tube Joint PTFE 테프론 튜브 조인트 마개

Cat. No.	Model	Connection For O.D	PK Price ₩
S12-460-011	D648-02	4 / 6 / (1/4") 6.35 mm	10ea 80,500
S12-460-015	D648-08	8 / (3/8") 9.52 / 10 mm	5ea 53,000
S12-460-019	D648-14	12 / (1/2") 12.7 / 14 / 16 mm	3ea 56,300

- 재질 : PTFE, -200 °C to 250°C 내열
- 테프론 재질로 최상의 내화학성을 가짐
- BOLA

