

Adapter, with Hose Connection (Cone Type)
상부 호스 연결 어댑터, LukeGL®

- For Tubing I.D 8mm
- Borosilicate Glass 3.3
- 요청에 의해 규격 외 사이즈도 공급 가능
- LukeGL®

일자형, Straight ①

Cat. No.	Model	Cone	EA Price ₩
A01-83-088	5092001	14/23	11,100
A01-83-096	5092003	24/40	16,000
A01-83-100	5092004	29/32	17,200

곡형, 90° Bent, Type I ②


Cat. No.	Model	Cone	EA Price ₩
A01-83-104	5095001	14/23	12,000
A01-83-112	5095003	24/40	16,000
A01-83-116	5095004	29/32	17,200

곡형, 90° Bent, Type II ③

Cat. No.	Model	Cone	EA Price ₩
A01-83-201	5095101	14/23	12,000
A01-83-205	5095103	24/40	16,000
A01-83-209	5095104	29/32	17,200

2-Side, 90° Bent x 2개 ④

Cat. No.	Model	Cone	EA Price ₩
A01-83-207	5095301	24/40	17,200



Adapter with Hose connection (Socket Type)
하부 호스 연결 어댑터, LukeGL®


- Borosilicate Glass 3.3
- LukeGL®

일자형, Straight ①

Cat. No.	Model	Socket	EA Price ₩
A01-83-120	5093001	14/23	11,100
A01-83-128	5093003	24/40	16,000
A01-83-132	5093004	29/32	17,200

곡형, Bent ②

Cat. No.	Model	Socket	EA Price ₩
A01-83-136	5096001	14/23	12,000
A01-83-144	5096003	24/40	16,000
A01-83-148	5096004	29/32	17,200



Adapter, with Hose Connection
호스 연결 어댑터

- O.D 10 mm Tubing과 같이 사용하기에 알맞음
- ASTM E438, Type I, Class A Glass
- Chemglass

일자형, Straight, Cone Type ①

Cat. No.	Model	Socket	EA Price ₩
A01-132-490	1012-14	14/20	42,300
A01-132-492	1012-19	19/22	44,200
A01-132-494	1012-01	24/40	47,000
A01-132-496	1012-02	29/42	53,600

곡형, 90° Bent, Cone Type ②

Cat. No.	Model	Cone	EA Price ₩
A01-132-421	1014-14	14/20	44,200
A01-132-423	1014-19	19/22	46,300
A01-132-427	1014-01	24/40	48,000
A01-132-429	1014-02	29/42	49,200

