

PC Container with Cover 투명 PC 원형 용기

Cat. No.	Model	용량	Size	EA Price ₩
C21-517-002	PCC02QT	1.8 L	Ø180×h105 mm	11,100
C21-517-004	PCC04QT	3.7 L	Ø180×h210 mm	11,600
C21-517-006	PCC06QT	5.6 L	Ø225×h195 mm	15,500
C21-517-008	PCC08QT	7.5 L	Ø225×h270 mm	19,800
C21-517-012	PCC12QT	11.3 L	Ø315×h190 mm	23,700
C21-517-018	PCC18QT	17 L	Ø315×h300 mm	41,500
C21-517-022	PCC22QT	20 L	Ø315×h380 mm	51,300

- 투명 Polycarbonate 재질
- Molded Graduation & Handle
- Rayhope®



Sample Cup for Tecnicon-Analyzer 샘플 컵

Cat. No.	Model	용량	Size (상부/하부 Ø × h)	PK Price ₩
C21-59-301	115015	1.5 ml	15/12.2×22.7 mm	1000ea 58,700
C21-59-303	115016	2 ml	14.8/12.8×24.9 mm	1000ea 60,800
C21-59-305	115017	4 ml	17/13.3×38 mm	1000ea 83,100

- 재질 : 투명 PS, -20 ~ 70°C 내열
- CE-marked according to IVD-Directive 98/79 EC
- Brand GmBH

Push-on Cap for #115015 & #115016

Cat. No.	Model	재질	적용 Cup	PK Price ₩
C21-59-309	115020	PE	for 1.5 / 2 ml Cup	1000ea 60,000



Sample Cup for Leuko-Ery Counter 샘플 컵 ①

Cat. No.	Model	용량	Size (Ø × h)	PK Price ₩	CS Price ₩
C21-59-323	722060	12 ml	22×38 mm	100ea 38,000	10pk 355,100

- 재질 : PS Body, PE Snap Cap
- Brand GmBH



Sample Cup for Coulter Counter 샘플 컵 ②

Cat. No.	Model	용량	Size (Ø × h)	PK Price ₩	CS Price ₩
C21-59-321	722055	20 ml	32×56 mm	100ea 26,600	10pk 248,300

- 재질 : PS Body, PE Lid
- CE-marked according to IVD-Directive 98/79 EC
- Brand GmBH

Planton-Plant Tissue Culture Container 식물 조직 배양 컨테이너

Cat. No.	Model	재질	EA Price ₩	PK Price ₩
C21-46-111	T20080	PC / PP	4,600	12ea 52,000

- 재질 : PS Body / PP Cover, Autoclavable
- Unique Cap and Drip Ring Shape
- Ideal for Plant Tissue Culture
- Size : w75×d75×h100 mm
- Tarsons

