

Hydrolyzing and Digestion Vessel for Microwave Oven

마이크로 오븐용 용기, 가수분해 다이제스트용

Cat. No.	Model	용량	Size	EA Price ₩
C21-460-011	A240-02	5 ml	∅15 × h32 mm	826,100
C21-460-019	A240-06	20 ml	∅22 × h60 mm	898,900
C21-460-023	A240-08	50 ml	∅33 × h62 mm	1,382,900
C21-460-027	A240-10	100 ml	∅35 × h110 mm	1,345,400

- 재질 : TFM , -200~250℃ 내열
 - Modified PTFE with Thermoplastic parts for Homogenous, Non-porous surface which avoids contamination and Memory effects
- 구성 : TFM Body and Screw Cover, One set of sealing and rupture membranes (10개)
- Quick and Easy Digestion용
- Microwave Oven의 Energy가 거의 손실없이 Vessel 안의 물질에 전달 됨
- 수조안에 Vessel 속의 용액을 Boiling Point 이상으로 가열 함
- Max. Pressure가 넘어가면 Membrane이 터지면서 Release된 Liquid를 Optional Tubing을 통해 Collection Vessel로 받을 수 있음
- BOLA

Model	O.D of Body	O.D of cover	Pressure Max. bar	Temperature Max. °C
A240-02	30 mm	40 mm	25 bar	160 °C
A240-06	50 mm	60 mm	20 bar	160 °C
A240-08	69 mm	95 mm	20 bar	160 °C
A240-10	70 mm	95 mm	15 bar	160 °C



A240-02

- PFA Vessel은 149 Page 참조
- PFA Sample Bottle은 150 Page 참조

