

Cylindrical Stainless Steel Container 타원형 스테인레스 용기

Cat. No.	Model	용량	Size	EA Price ₩
C21-177-010	LCC-10	1.1 L	Ø100×h150 mm	14,300
C21-177-020	LCC-12	1.7 L	Ø120×h150 mm	17,000
C21-177-030	LCC-15	3.1 L	Ø150×h180 mm	20,500
C21-177-040	LCC-16	3.8 L	Ø160×h190 mm	23,400
C21-177-050	LCC-18	5.3 L	Ø180×h210 mm	28,100
C21-177-060	LCC-20	7.5 L	Ø200×h240 mm	33,200
C21-177-070	LCC-23	13.2 L	Ø230×h270 mm	39,700

- LK Labkorea



Stainless Steel Container with Hole 타공 스테인레스 용기

Cat. No.	Model	용량	Size	EA Price ₩
C21-177-361	LCH-055	7 L	Ø250×h150 mm	67,500
C21-177-362	LCH-056	12 L	Ø250×h250 mm	79,600

- 재질 : Stainless Steel 304, Silicon Gasket
- Clip 및 Cover 포함
- LK Labkorea



Square Stainless Steel Container 사각 스테인레스 용기

Cat. No.	Model	용량	Size (w×d×h, mm)	EA Price ₩
C21-177-391	SSC-001	4.5 L	230×170×120	53,100
C21-177-392	SSC-002	6.5 L	260×200×130	62,700
C21-177-393	SSC-003	9.5 L	290×230×150	72,300
C21-177-394	SSC-004	12.5 L	310×260×160	84,400

- 재질 : Stainless Steel 304, Silicon Gasket
- Clip 및 Cover 포함
- LK Labkorea



Double Walled Stainless Steel Container 이중 스테인레스 용기

Cat. No.	Model	용량	Size	EA Price ₩
C21-177-371	DWC-001	10 L	Ø290×h180 mm	127,800
C21-177-373	DWC-003	16 L	Ø290×h280 mm	144,600
C21-177-374	DWC-004	18 L	Ø290×h320 mm	156,700
C21-177-375	DWC-005	12 L	Ø325×h180 mm	127,800
C21-177-377	DWC-007	18 L	Ø325×h280 mm	156,700
C21-177-379	DWC-009	21 L	Ø325×h320 mm	168,700
C21-177-380	DWC-010	30 L	Ø370×h350 mm	216,900

- 재질 : Stainless Steel 304, Silicon Gasket
- 이중 타입의 스테인레스 원형 용기
- LK Labkorea

