

## Stainless Steel Container 스테인레스 용기

- 재질 : Stainless Steel 304, Silicon Gasket
- 용기와 뚜껑의 결합 방식으로 구분   • Buff Polish
- LK Labkorea

### Standard Type

Cat. No.	Model	용량	Size	배송비 별도	
				EA Price	₩
C21-177-111	LSCG-10	10 L	Ø250×h250 mm	72,300	
C21-177-115	LSCG-20	20 L	Ø300×h300 mm	96,400	
C21-177-119	LSCG-40	40 L	Ø350×h390 mm	144,600	
C21-177-123	LSCG-60	60 L	Ø400×h460 mm	192,800	
C21-177-127	LSCG-80	80 L	Ø450×h480 mm	209,700	
C21-177-131	LSCG-100	100 L	Ø490×h510 mm	241,000	
C21-177-135	LSCG-120	120 L	Ø550×h510 mm	265,100	
C21-177-139	LSCG-160	160 L	Ø550×h630 mm	289,200	
C21-177-143	LSCG-200	200 L	Ø635×h640 mm	361,500	
C21-177-147	LSCG-260	260 L	Ø635×h810 mm	409,700	

- 구성 : 본체 + 뚜껑

### Clip Type

Cat. No.	Model	용량	Size	배송비 별도	
				EA Price	₩
C21-177-211	LSCC-10	10 L	Ø250×h250 mm	108,500	
C21-177-215	LSCC-20	20 L	Ø300×h300 mm	120,500	
C21-177-219	LSCC-40	40 L	Ø350×h390 mm	178,200	
C21-177-223	LSCC-60	60 L	Ø400×h460 mm	204,900	
C21-177-227	LSCC-80	80 L	Ø450×h480 mm	229,000	
C21-177-231	LSCC-100	100 L	Ø490×h510 mm	277,200	
C21-177-235	LSCC-120	120 L	Ø550×h510 mm	301,300	
C21-177-239	LSCC-160	160 L	Ø550×h630 mm	325,400	
C21-177-243	LSCC-200	200 L	Ø635×h640 mm	482,000	
C21-177-247	LSCC-260	260 L	Ø635×h810 mm	753,100	

- 구성 : 본체 + 뚜껑 + 클립 + Silicon Gasket

### Clamp Type

Cat. No.	Model	용량	Size	배송비 별도	
				EA Price	₩
C21-177-311	LSCL-10	10 L	Ø250×h250 mm	216,900	
C21-177-315	LSCL-20	20 L	Ø300×h300 mm	229,000	
C21-177-319	LSCL-30	30 L	Ø350×h390 mm	265,100	
C21-177-323	LSCL-50	50 L	Ø400×h460 mm	337,400	
C21-177-327	LSCL-80	80 L	Ø450×h480 mm	337,400	
C21-177-331	LSCL-100	100 L	Ø490×h510 mm	385,600	
C21-177-335	LSCL-120	120 L	Ø550×h510 mm	409,700	
C21-177-339	LSCL-160	160 L	Ø550×h630 mm	433,800	
C21-177-343	LSCL-200	200 L	Ø635×h640 mm	506,100	
C21-177-347	LSCL-260	260 L	Ø635×h810 mm	554,300	

- 구성 : 본체 + 뚜껑 + 클램프 + Silicon Gasket
- 밀폐율이 가장 우수함



## Tight Sealing Stainless Steel Container 밀폐형 스테인레스 용기

Cat. No.	Model	용량	Size	Thickness	EA Price	₩
C21-215-010	5-152-06	10 L	Ø240×h240 mm	0.7 mm	819,700	
C21-215-020	5-152-05	20 L	Ø300×h300 mm	0.8 mm	1,021,300	

- 재질 : Stainless Steel #304, Silicon Rubber, PC Window
- 뚜껑에 투명한 창이 있어서 열지 않고 내부 확인 가능
- Buff Polish
- Universal

