

Digital Colony Counter 디지털 세균 계수기

Cat. No.	Model	Description	EA Price ₩
C23-108-022	3081502	Colony Counter Set	3,335,200

- 배양된 대장균 및 일반 세균등의 계수용 • “9999” 까지 기록할 수 있음
- 본체 위, 아래, 사이드 부분으로 부터의 3가지 조명 모드 가능 (위는 Option임)
- 평균값 계산 (최대 99회, 32,300 Colonies)
- 세균에 열 변성을 가져오지 않는 LED Light를 통해 Colony의 높은 투명도를 보장
- Sound와 Count가 동시에 이루어짐 • USB를 통해 데이터 전송가능
- 구성 : 본체, 90 mm용 Petridish Holder, Marking Pen Probe,
Pen Holder Contrast Disk (Black)
- Schuett


Cat. No.	C23-108-022
Model	3081502
Display	4 Digital LED Display
Petri Dish Holder	90~100 mm용 기본, 50-60 mm와 140-150 mm은 Option
Pressure Sensitivity	Adjustable
Size	w280×d130×h250 mm
Mode of Counting	Needle - Point Probing, Marking Pen Probing
Power	100 ~ 240 V, 50/60 Hz

Optional Accessory

Cat. No.	Model	Description	EA Price ₩
C23-108-522	3081702	LED top Light Lamp, Flexible	514,100
C23-108-546	3081802	Adapter for Petridishes with 50-60mm	123,500
C23-108-569	3081812	Adapter for Petridishes with 140-150mm	227,100
C23-108-575	3081922	Counting Needle	96,700
C23-108-589	3081932	Felt-tip Marker Pen	15,100
C23-108-610	3081602	3×Magnifying Glass, Flexible	313,900
C23-108-625	3081902	Spiral-Plater Disk, Black/White	140,100



Illumination




Automatic Colony Counter 자동 세균 계수기


Cat. No.	Model	Description	EA Price ₩
C23-108-550	ColonyQuant	Automatic Colony Counter Set	

- 콜로니의 Color, Size and Shape에 따라 Software를 선택 가능함
- Light - Proof Sample Chamber with CCD - Firewire - Camera (색상)
- Image Acquisition and Simultaneous Evaluation within Seconds
- 테이블과 이미지가 디지털로 저장됨
- Agar Plates, Nutrient Disks, Inhibition Zone Analysis, Spiral Plating
- Mixed Cultures, up to 8 Colours Simultaneously per Petridish
- Sectorial Evaluation, 90% Option
- In Conformity with GLP • 기타 자세한 사항은 문의 요함
- Schuett


Cat. No.	C23-108-550
Model	Colony Quant
Resolution	640×480 pixel
Illumination	From Below or Above Selectable
Size	w240×d290×h479 mm
Feature	Light - Proof Sample Chamber
Power	100 ~ 240 V, 50/60 Hz



Enterobacteria




Brewer's yeast



E.coli colonies growing on LB-agar



Total numbers of colonies



E.coli colonies