

Alumina Crucible 알루미나 도가니

- 알루미나(Al_2O_3) 99.8%
- 내열온도 : 1,700°C
- 최상의 내화학성, 내열성, 내마모성
- LabCeramic®

Cat. No.	Model	용량	Size	EA Price ₩
C24-45-714	LAC1926S	5 ml	∅19×h26 mm	13,500
C24-45-721	LAC2929S	10 ml	∅29×h30 mm	14,300
C24-45-735	LAC4030S	20 ml	∅39×h29 mm	13,600
C24-45-742	LAC4637S	30 ml	∅52×h39 mm	17,700
C24-45-756	LAC5745S	50 ml	∅61×h49 mm	26,300
C24-45-763	LAC7055S	100 ml	∅72×h58 mm	36,300
C24-45-770	LAC7658S	150 ml	∅81×h67 mm	43,900
C24-45-777	LAC8164S	200 ml	∅86×h70 mm	47,300
C24-45-784	LAC8574S	250 ml	∅92×h75 mm	52,300
C24-45-791	LAC9280S	300 ml	∅97×h76 mm	64,000
C24-45-798	LAC1110S	500 ml	∅110×h105 mm	109,600
C24-45-805	LAC1412S	1 L	∅140×h117 mm	195,700



Body Only

Cat. No.	Model	용량	EA Price ₩
C24-45-014	LAC1926	5 ml	7,800
C24-45-021	LAC2929	10 ml	8,800
C24-45-035	LAC4030	20 ml	7,100
C24-45-042	LAC4637	30 ml	10,600
C24-45-056	LAC5745	50 ml	17,600
C24-45-063	LAC7055	100 ml	26,400
C24-45-070	LAC7658	150 ml	33,400
C24-45-077	LAC8164	200 ml	35,200
C24-45-084	LAC8574	250 ml	38,600
C24-45-091	LAC9280	300 ml	47,300
C24-45-098	LAC1110	500 ml	91,200
C24-45-105	LAC1412	1 L	166,700

Cover Only

Cat. No.	Model	Description	EA Price ₩
C24-45-514	LAC1926C	5 ml 용	6,100
C24-45-521	LAC2929C	10 ml 용	6,100
C24-45-535	LAC4030C	20 ml 용	7,100
C24-45-542	LAC4637C	30 ml 용	7,800
C24-45-556	LAC5745C	50 ml 용	9,800
C24-45-563	LAC7055C	100 ml 용	11,400
C24-45-570	LAC7658C	150 ml 용	12,300
C24-45-577	LAC8164C	200 ml 용	14,100
C24-45-584	LAC8574C	250 ml 용	15,900
C24-45-591	LAC9280C	300 ml 용	19,300
C24-45-598	LAC1110C	500 ml 용	22,900
C24-45-605	LAC1412C	1 L 용	36,900

Alumina Cylindrical Crucible 알루미나 원통형 도가니

- 알루미나(Al_2O_3) 99.8%
- 내열온도 : 1,700°C
- 최상의 내화학성, 내열성, 내마모성
- LabCeramic®

Body with Cover Set

Cat. No.	Model	용량	Size	EA Price ₩
C24-45-815	ACC001S	1 ml	∅13×h21.5 mm	20,200
C24-45-820	ACC002S	2 ml	∅13×h31.5 mm	21,800
C24-45-825	ACC028S	28 ml	∅36×h48.5 mm	14,400
C24-45-830	ACC060S	60 ml	∅44×h63.5 mm	27,000
C24-45-835	ACC170S	170 ml	∅56×h98.5 mm	48,400
C24-45-840	ACC340S	340 ml	∅72×h108.5 mm	75,200
C24-45-845	ACC400S	400 ml	∅76×h113.5 mm	81,700



Body Only

Cat. No.	Model	용량	Size	EA Price ₩
C24-45-205	ACC001B	1 ml	∅13×h20 mm	10,600
C24-45-210	ACC002B	2 ml	∅13×h30 mm	12,300
C24-45-215	ACC028B	28 ml	∅36×h45 mm	10,200
C24-45-220	ACC060B	60 ml	∅44×h60 mm	19,300
C24-45-225	ACC170B	170 ml	∅56×h95 mm	35,200
C24-45-230	ACC340B	340 ml	∅72×h105 mm	55,700
C24-45-235	ACC400B	400 ml	∅76×h110 mm	75,900

Cover Only

Cat. No.	Model	용량	Size	EA Price ₩
C24-45-905	ACC012C	1 ml, 2 ml	∅13 mm	10,600
C24-45-915	ACC028C	28 ml	∅36 mm	5,400
C24-45-920	ACC060C	60 ml	∅44 mm	8,800
C24-45-925	ACC170C	170 ml	∅56 mm	18,800
C24-45-930	ACC340C	340 ml	∅72 mm	23,500
C24-45-935	ACC400C	400 ml	∅76 mm	25,100