


Automatic Film Applicator 자동 도포기

Cat. No.	Model	Bed Type	Bed Size (mm)	EA Price ₩
A03-178-010	AFA350	Glass	w300×d490	3,578,400
A03-178-020	AFA350V	Vacuum		4,879,600
A03-178-030	AFA350VH	Vacuum & Heating		7,048,200


- 페인트, 잉크, 수지등 정확 재현성 있는 샘플 코팅
- 재현 가능한 모든 요소들 제공 (속도, 방향, 압력)
- 수동 코팅시 생기는 인적 오류 없이 항상 일정한 결과치를 얻을 수 있음
- 주요기능 : 시험속도, 이동 거리 조절 가능
- Heating Type : Max.150 °C
- 잔고장 없는 견고한 디자인
- 용도에 따라 Wire Bar Coater 또는 Micro Applicator와 같이 사용 가능함
- 용도에 따라 필요한 Bed로 교체 가능함
- Vacuum Bed에 사용되는 Pump는 옵션임
- LK Labkorea



Cat. No.	A03-178-010	A03-178-020	A03-178-030
Model	AFA350	AFA350V	AFA350VH
Bed Type	Glass	Vacuum	Vacuum & Heating
Plate 재질	Glass	Aluminum Anodizing	
Size	Bed	w300×d490 mm	
	Coating	w270×d380 mm	
	Overall	w680×d500×h360 mm	
Speed	10 ~ 100 mm/sec (Digital Display)		
Controller	Digital Multi-Function Controller		
Safety	Limit Switch & 위치 설정 센서		
Power	220V, 60Hz		

Wire Bar Coater 바코타, “AFA350 Auto Flim Applicator” 용

Rod, Ø1/4" (6.35 mm)			Rod, Ø1/2" (12.7 mm)			Wet Film (µm)
Cat. No.	Model	EA Price ₩	Cat. No.	Model	EA Price ₩	
A03-43-008	LWBO003	93,700	A03-43-272	LWBO103	107,100	6.86
A03-43-016	LWBO004	93,700	A03-43-280	LWBO104	107,100	9.14
A03-43-024	LWBO005	93,700	A03-43-288	LWBO105	107,100	11.43
A03-43-032	LWBO006	93,700	A03-43-296	LWBO106	107,100	13.7
A03-43-040	LWBO007	93,700	A03-43-304	LWBO107	107,100	16.0
A03-43-048	LWBO008	93,700	A03-43-312	LWBO108	107,100	18.3
A03-43-056	LWBO009	93,700	A03-43-320	LWBO109	107,100	20.6
A03-43-064	LWBO010	93,700	A03-43-328	LWBO110	107,100	22.9
A03-43-072	LWBO012	93,700	A03-43-336	LWBO112	107,100	27.4
A03-43-080	LWBO014	93,700	A03-43-344	LWBO114	107,100	32.0
A03-43-088	LWBO016	93,700	A03-43-352	LWBO116	107,100	36.6
A03-43-096	LWBO018	93,700	A03-43-360	LWBO118	107,100	41.1
A03-43-104	LWBO020	93,700	A03-43-368	LWBO120	107,100	45.7
A03-43-112	LWBO022	93,700	A03-43-376	LWBO122	107,100	50.3
A03-43-120	LWBO024	93,700	A03-43-384	LWBO124	107,100	54.8
A03-43-128	LWBO026	93,700	A03-43-392	LWBO126	107,100	59.4
A03-43-136	LWBO028	93,700	A03-43-400	LWBO128	107,100	64.8
A03-43-144	LWBO030	93,700	A03-43-408	LWBO130	107,100	68.6
A03-43-152	LWBO032	93,700	A03-43-416	LWBO132	107,100	73.2
A03-43-160	LWBO034	93,700	A03-43-424	LWBO134	107,100	77.7
A03-43-168	LWBO036	93,700	A03-43-432	LWBO136	107,100	83.3
A03-43-176	LWBO038	93,700	A03-43-440	LWBO138	107,100	86.9
A03-43-184	LWBO040	93,700	A03-43-448	LWBO140	107,100	91.4
A03-43-192	LWBO042	93,700	A03-43-456	LWBO142	107,100	96.0
A03-43-200	LWBO044	93,700	A03-43-464	LWBO144	107,100	100.6
A03-43-208	LWBO046	93,700	A03-43-472	LWBO146	107,100	105.2
A03-43-216	LWBO048	93,700	A03-43-480	LWBO148	107,100	109.8
A03-43-224	LWBO050	93,700	A03-43-488	LWBO150	107,100	114.3
A03-43-232	LWBO055	93,700	A03-43-496	LWBO155	107,100	125.7
A03-43-240	LWBO060	93,700	A03-43-504	LWBO160	107,100	137.2
A03-43-248	LWBO065	93,700	A03-43-512	LWBO165	107,100	148.6
A03-43-256	LWBO070	93,700	A03-43-520	LWBO170	107,100	160.0
A03-43-264	LWBO075	93,700	A03-43-528	LWBO175	107,100	171.5



- Rod 지름 : Ø1/4", Ø1/2"
- Coating Width : 300 mm
- Overall 길이 : 400 mm
- 재질 : Stainless Steel 304
- LK Labkorea