

## Glass Filter Crucible, High Quality 고품질 도가니형 글라스 필터, ISO

Cat. No.	Model	용량	Porosity	Porosity 필터	Size	EA Price ₩
F07-228-216	20301	30 ml	1	Ø 30 mm	36×56 mm	22,400
F07-228-222	20302		2			22,400
F07-228-228	20303		3			22,400
F07-228-234	20304		4			25,800
F07-228-240	20305		5			47,300
F07-228-246	20501	50 ml	1	Ø 40 mm	46×61 mm	32,400
F07-228-252	20502		2			32,400
F07-228-258	20503		3			32,400
F07-228-264	20504		4			37,100
F07-228-270	20505		5			70,500

- Borosilicate Glass 3.3
- ISO 규격에 따른 Porosity
- Heavy Wall Type으로 VitraPOR® 고품질 Filter 사용
- Robu / VitraPOR®



Porosity	ISO4793	Pore Size
1	P160	100 ~ 600 um
2	P100	40 ~ 100 um
3	P40	16 ~ 40 um
4	P16	10 ~ 16 um
5	P1.6	1.0 ~ 1.6 um

P1	Coarse Filtration, Gas Distribution, Liquid Distribution, Gas Cleaning
P2	Preparative fine Filtration, Mercury Filtration, Preparative Work with Coarse Crystalline Products
P3	Analytical Work with Medium Fine Precipitates Analytical Filtration, Fine Gas Filtration, Special Filtration
P4	Analytical Fine Filtration, Analytical and Preparative Works with The Finest Precipitates Return and Stop Valves
P5	General Bacteria Filtration, Sterile Filtration

## Glass Filter Crucible 도가니형 글라스 필터, LukeGL®

Cat. No.	Model	용량	Porosity	Porosity 필터	EA Price ₩
C24-503-004	GFC130	30 ml	1	Ø 21 mm	12,300
C24-503-007	GFC230		2		12,300
C24-503-010	GFC330		3		12,300
C24-503-013	GFC430		4		12,300
C24-503-016	GFC530		5		12,300
C24-503-022	GFC150	50 ml	1	Ø 27 mm	19,600
C24-503-025	GFC250		2		19,600
C24-503-028	GFC350		3		19,600
C24-503-031	GFC450		4		19,600
C24-503-034	GFC550		5		19,600

- Borosilicate Glass 3.3
- LukeGL®



Porosity	Pore Size
1	90 ~ 150 um
2	40 ~ 90 um
3	15 ~ 40 um
4	5 ~ 15 um
5	1.0 ~ 2.0 um

## Filter Crucible Adapter 도가니형 글라스 필터 어댑터, LukeGL®

Cat. No.	Model	적용 Crucible	id	Height	EA Price ₩
C24-52-511	L369002M	10 ml	30 mm	60 mm	21,700
C24-52-515	L369004M	30 ml	37 mm	68 mm	22,500
C24-52-519	L369008M	60 ml	47 mm	78 mm	23,600

- Borosilicate Glass 3.3
- Stem Size : Ø9×L75 mm
- LukeGL®



- Filter Adapter는 493 Page 참조

