

## Porcelain Crucible

### 자제 도가니

- Max. Temp : 1,000°C
- Water Absorption 용량 : 0%
- Density : 2.4 g/cm<sup>3</sup>
- Jipo / LabCeramic®

#### Low Form / 단형 자제 도가니 Only ①

Cat. No.	Model	용량	Size	EA Price ₩
C24-128-111	1/30	5 ml	∅30×h19 mm	3,500
C24-128-113	1/35	10 ml	∅35×h22 mm	4,000
C24-128-115	1/40	17 ml	∅40×h25 mm	4,100
C24-128-117	1/45	21 ml	∅45×h28 mm	5,200
C24-128-119	1/50	34 ml	∅50×h32 mm	5,700
C24-128-121	1/60	62 ml	∅60×h38 mm	7,800
C24-128-123	1/70	91 ml	∅70×h44 mm	10,500

#### Middle Form / 중형 자제 도가니 Only ②

Cat. No.	Model	용량	Size	EA Price ₩
C24-128-131	2/30	10 ml	∅30×h25 mm	3,500
C24-128-133	2/35	12 ml	∅35×h28 mm	3,400
C24-128-135	2/40	20 ml	∅40×h32 mm	4,100
C24-128-137	2/45	30 ml	∅45×h36 mm	4,100
C24-128-139	2/50	45 ml	∅50×h40 mm	5,400
C24-128-141	2/60	80 ml	∅60×h48 mm	6,700
C24-128-143	2/70	120 ml	∅70×h56 mm	9,000
C24-128-145	2/80	200 ml	∅80×h62 mm	9,600

#### High Form / 장형 자제 도가니 Only ③

Cat. No.	Model	용량	Size	EA Price ₩
C24-128-161	3/30	15 ml	∅30×h38 mm	5,200
C24-128-163	3/35	26 ml	∅35×h44 mm	6,400
C24-128-165	3/40	35 ml	∅40×h50 mm	7,100
C24-128-167	3/45	50 ml	∅45×h56 mm	7,500
C24-128-169	3/50	72 ml	∅50×h62 mm	8,000
C24-128-171	3/60	120 ml	∅60×h75 mm	9,400

#### Crucible Lid / 자제 도가니 뚜껑 Only ④

Cat. No.	Model	Size	사용 도가니	EA Price ₩
C24-128-201	D30	od34 mm	od30 mm용	3,100
C24-128-203	D35	od39 mm	od35 mm용	3,400
C24-128-205	D40	od44 mm	od40 mm용	3,600
C24-128-207	D45	od49 mm	od45 mm용	4,500
C24-128-209	D50	od54 mm	od50 mm용	4,800
C24-128-211	D60	od64 mm	od60 mm용	5,300
C24-128-213	D70	od74 mm	od70 mm용	6,600
C24-128-215	D80	od84 mm	od80 mm용	7,000

## Porcelain Gooch Filtering Crucible

### 자제 구찌 도가니

Cat. No.	Model	용량	Size	EA Price ₩
C24-128-313	234/1	17 ml	∅30×h36 mm	4,800
C24-128-315	234/3	25 ml	∅35×h43 mm	5,200
C24-128-317	234/4	35 ml	∅39×h49 mm	5,800
C24-128-319	234/5	120 ml	∅60×h71 mm	7,100

- Max. Temp : 1,000°C
- Water Absorption 용량 : 0 %
- Density : 2.4 g/cm<sup>3</sup>
- Jipo / LabCeramic®

