

Nickel Crucible with Cover

니켈 도가니

Cat. No.	Model	용량	Size	EA Price ₩
C24-102-015	L810484	30 ml	∅38×h38 mm	43,900
C24-102-020	L810585	50 ml	∅48×h43 mm	50,100

- Nickel 99.5%
- 내열 온도 : 800℃
- LK Labkorea



Platinum Crucible with Cover

백금 도가니

Cat. No.	Model	용량	Size	EA Price ₩
C24-103-008	PC0010S	10 ml	∅25×h27 mm	
C24-103-016	PC0015S	15 ml	∅28×h31 mm	
C24-103-024	PC0020S	20 ml	∅32×h34 mm	
C24-103-032	PC0025S	25 ml	∅34×h38 mm	
C24-103-040	PC0030S	30 ml	∅36×h40 mm	
C24-103-048	PC0040S	40 ml	∅41×h45 mm	
C24-103-056	PC0050S	50 ml	∅44×h47 mm	
C24-103-064	PC0070S	70 ml	∅48×h51 mm	
C24-103-072	PC0100S	100 ml	∅56×h62 mm	
C24-103-080	PC0150S	150 ml	∅62×h68 mm	
C24-103-088	PC0200S	200 ml	∅68×h74 mm	
C24-103-096	PC0250S	250 ml	∅72×h82 mm	

- Platinum 99.95%
- 내열 온도 : 1,700℃, 적정 온도 : 800 ~ 1,400℃
- 기본적으로 뚜껑이 포함됨
- 요청에 의해 뚜껑없이 같은 무게로 도가니 하부로만 공급 가능
- 백금 도가니 용량 = 백금 무게 (시세에 따라 가격 변동)
- LK Labkorea




PTFE Crucible

PTFE 테프론 도가니

Cat. No.	Model	용량	Size	EA Price ₩
D02-59-222	009.005	5 ml	∅27×h16 mm	14,200
D02-59-226	009.025	25 ml	∅44×h26 mm	22,000
D02-59-230	009.075	75 ml	∅56×h41 mm	28,600
D02-59-234	009.100	100ml	∅60×h45 mm	40,700

- 재질 : PTFE, ~200 ~ 280℃ 내열
- 최상의 내화학성, 내열성, 내마모성
- Crucible Form Without Spout
- Cowie



Quartz Crucible

석영 도가니

Cat. No.	Model	용량	Size	EA Price ₩
C24-60-050	QC0015	15 ml	∅38, ∅25×h28 mm	30,000
C24-60-150	QC0030	30 ml	∅47, ∅25×h36 mm	40,300
C24-60-200	QC0050	50 ml	∅54, ∅34×h42 mm	56,000
C24-60-250	QC0100	100 ml	∅70, ∅34×h55 mm	82,000

- 내열 온도 : 1,250℃, 상용 온도 : 1,100℃
- 뚜껑 포함
- LK Labkorea

