

Stainless Steel Crucible 스테인레스 도가니, without Cover

Cat. No.	Model	용량	Size	EA Price ₩
C24-284-050	5640	35 ml	∅40×h40 mm	19,900
C24-284-100	5641	45 ml	∅45×h40 mm	21,800
C24-284-150	5642	65 ml	∅50×h45 mm	23,600
C24-284-200	5643	100 ml	∅60×h50 mm	25,500

- 재질 : Stainless Steel 18/8 , Thickness : 1 mm
- Usbeck



Silicon Carbide Crucible 실리콘 카바이드 도가니


- 재질 : 99% Sic, 0.2% SiO₂, 0.1% C
- Max. Temp : 1600°C
- 최상의 내화학성, 내열성, 내마모성
- Cover는 별도 구매
- Universal

Body Only

Cat. No.	Model	용량	Size	EA Price ₩
C24-215-030	1-9669-01	100 ml	∅60×h70 mm	493,000
C24-215-032	1-9669-02	500 ml	∅100×h115 mm	964,500

Cover Only

Cat. No.	Model	용량	Size	EA Price ₩
C24-215-034	1-9669-11	100 ml	∅63×h6 mm	321,500
C24-215-036	1-9669-12	500 ml	∅103×h6 mm	460,800



Zirconia Crucible SET 지르코니아 도가니 세트

Cat. No.	Model	용량	Size	EA Price ₩
C24-45-941	LZC039029S	20ml	∅39 x H33 mm	82,700
C24-45-943	LZC049036S	30ml	∅49 x H39 mm	154,300
C24-45-945	LZC057044S	50ml	∅57 x H47 mm	236,100
C24-45-947	LZC072059S	100ml	∅72 x H62 mm	283,500
C24-45-949	LZC081065S	150ml	∅81 x H68 mm	370,200
C24-45-951	LZC090069S	250ml	∅90 x H73 mm	750,600
C24-45-953	LZC110101S	500ml	∅110 x H105 mm	846,700

- 산화지르코늄(ZrO₂) 94.2% • 내열 온도 : 1700°C
- 1°C/min 이하의 가열 및 냉각조건 필수 • 고온, 내화학 환경에 적합
- LabCeramic®



- 사각, 원형 알루미늄 Dish는 348 Page 참조

- 연소용 자재 Dish는 348 Page 참조