

## Wide Mouth Dewar Flask

### 광구형 드와 플라스크

- Vented PE Stopper 포함
- 가장 기본적인 Dewar Flask
- 외부 열전달 차단을 위해 Dewar Flask 사이에 패드가 없음
- 메쉬 커버 사용 : 그립감 및 안정성 증대
- Without Handle : Protective Plastic Cover 포함
- POPE

#### Without Handle ①

Cat. No.	Model	용량	Size (mm)		EA Price ₩
			Inside	Overall	
C25-411-131	8600	350 ml	Ø70×125	Ø86×157	259,300
C25-411-135	8640	665 ml	Ø70×195	Ø86×241	274,500
C25-411-140	8645	1 L	Ø70×302	Ø86×337	289,800
C25-411-145	8621	1.9 L	Ø120×190	Ø150×252	731,900
C25-411-150	8642	4.3 L	Ø143×305	Ø178×365	838,700

- Case 재질 : Metal Housing & Weighted Metal Base
- Protective Plastic Cover 포함

#### With Handle ②

Cat. No.	Model	용량	Size (mm)		EA Price ₩
			Inside	Overall	
C25-411-155	8945	1 L	Ø70×302	Ø86×337	396,500
C25-411-160	8921	1.9 L	Ø120×190	Ø150×252	854,000
C25-411-165	8942	4.3 L	Ø143×305	Ø178×365	945,400

- Case 재질 : Aluminum

#### Spare Stopper

Cat. No.	Model	Description	EA Price ₩
C25-411-901	8821	8621, 8921용 Stopper	22,900
C25-411-904	8842	8642, 8942용 Stopper	22,900
C25-411-907	8845	8600, 8640, 8645, 8945용 Stopper	16,700



## Spherical Dewar Flask

### 세구형 드와 플라스크

Cat. No.	Model	용량	Size (mm)		EA Price ₩
			Inside	Overall	
C25-411-180	8697	2 L	Ø158×335	Ø184×386	3,303,900
C25-411-185	8698	5 L	Ø216×370	Ø255×456	3,558,000

- Case 재질 : Aluminum
- Vented PE Stopper 포함
- 액체질소의 임시 저장 및 이송용
- 외부 열전달 차단을 위해 Dewar Flask 사이에 패드가 없음
- POPE

#### Narrow Mouth Caddle Only

Cat. No.	Model	Description	EA Price ₩
C25-411-190	30270	8697 용, with Handle	548,900
C25-411-195	30380	8698 용, with Handle	671,000

- Spherical Dewar Flask용 보관 및 이동 Case

#### Spare Stopper

Cat. No.	Model	Description	EA Price ₩
C25-411-200	7CORK-97	8657 용, Cork 재질	13,500
C25-411-202	7CORK-98	8698 용, Cork 재질	15,900

