



Liquid Nitrogen Storage Vessels

액체 질소 보관용 베셀

배송비 별도

Cat. No.	Model	용량	Size	EA Price ₩
C25-07-005	Cryolab 2	2 L	Ø185×h385 mm	2,009,000
C25-07-010	Cryolab 5	5 L	Ø216×h495 mm	2,671,000
C25-07-015	Cryolab 10	10 L	Ø305×h560 mm	2,704,500
C25-07-020	Cryolab 20	20 L	Ø407×h570 mm	3,174,800
C25-07-025	Cryolab 25	25 L	Ø395×h585 mm	3,158,100
C25-07-030	Cryolab 35	35 L	Ø407×h650 mm	3,346,100
C25-07-035	Cryolab 50	50 L	Ø407×h810 mm	4,014,800

- Strong, Lightweight Aluminium Construction
- 높은 열효율과 매우 낮은 증발 손실
- Unique Pouring Trolley for Cryolab 25
- 버킷 핸들 (2L, 5L & 10L) 또는 이중 핸들 (20L, 25L, 35L & 50L)
- Roller Base, Trolley 와 LN2 Dispenser는 별도 (Accessory란 참조)
- Statebourne / Cryolab®



Cat. No.	C25-07-005	C25-07-010	C25-07-015	C25-07-020	C25-07-025	C25-07-030	C25-07-035
Model	Cryolab 2	Cryolab 5	Cryolab 10	Cryolab 20	Cryolab 25	Cryolab 35	Cryolab 50
Capacity	2 L	5 L	10 L	20 L	25 L	35 L	50 L
Evaporation Rate	0.08 L/day	0.18 L/day	0.2 L/day	0.18 L/day	0.25 L/day	0.27 L/day	0.28 L/day
Static Hoding Time	25 day	28 day	50 day	111 day	100 day	130 day	178 day
Weight Empty	2.8 kg	4.1 kg	5.6 kg	9.1 kg	7.4 kg	12.2 kg	15.5 kg
Weight Full	4.4 kg	8.1 kg	13.7 kg	25.3 kg	27.6 kg	40.5 kg	55.9 kg
Neck Diameter	32 mm	51 mm	51 mm	51 mm	51 mm	63 mm	63 mm
Roller Base	×	×	×	○	○	○	○
LN ₂ Dispenser	×	×	○	○	○	○	○
LN ₂ Dispenser w/ Flexible Hose	×	×	×	○	○	○	○

Liquid Dispenser / Liquid Withdrawal Device ①

Cat. No.	Model	적용 Vessel Neck Size	EA Price ₩
C25-07-610	9701306	51 mm for Cryolab 10	2,106,300
C25-07-612	9701056	51 mm for Cryolab 20, 25	2,285,800
C25-07-614	9701057A	63 mm for Cryolab 35	2,362,900
C25-07-616	9701057	63 mm for Cryolab 50	2,362,900

* Liquid Withdrawal Device

Liquid Dispenser with 1.2m Flexible Hose ②

Cat. No.	Model	적용 Vessel Neck Size	EA Price ₩
C25-07-622	9701157	51 mm for Cryolab 20, 25 (w/ 1.2m Transfer Hose)	2,553,000
C25-07-624	9701152A	63 mm for Cryolab 35 (w/ 1.2m Transfer Hose)	2,629,900
C25-07-628	9701152	63 mm for Cryolab 50 (w/ 1.2m Transfer Hose)	2,629,900

* Liquid Withdrawal Device with 1.2m Flexible Hose

Roller Base ③

Cat. No.	Model	Description	EA Price ₩
C25-07-500	9902082	for Cryolab 20, 35, 50	375,700
C25-07-510	9901063	for Cryolab 25	375,700

Trolley for Cryolab 25 ④

Cat. No.	Model	Description	EA Price ₩
C25-07-515	9901064	Pouring Trolley, w472 x h1033 mm	1,186,000

