

Rack-type Cryosystem 랙 타입 액체 질소 저장 탱크

배송비 별도

Cat. No.	Model	용량	Size	EA Price ₩
C25-07-080	BioRack 750	35 L	Ø480×h747 mm	8,795,200
C25-07-085	BioRack 3000	85 L	Ø646×h738 mm	13,207,300
C25-07-087	BioRack 4800	131 L	Ø646×h913 mm	16,520,200
C25-07-090	BioRack 6000	160 L	Ø648×h1027 mm	17,811,900

- Strong, Lightweight Aluminium Construction
- High Thermal Efficiency
- 2 ml용 Racked Storage System
- Rack과 Cryobox가 기본 포함
- Long Term Storage를 위한 Autofill Accessory Option은 문의 요망
- Statebourne / BioRack®



Cat. No.	C25-07-080	C25-07-085	C25-07-087	C25-07-090
Model	BioRack 750	BioRack 3000	BioRack 4800	BioRack 6000
용량	35 L	85 L	131 L	160 L
Evaporation Rate	0.30 L/day	0.63 L/day	0.63 L/day	0.63 L/day
Static Holding Time	117 day	135 day	208 day	254 day
Weight Empty	19 kg	42 kg	46 kg	64 kg
Weight Full	47 kg	110 kg	152 kg	193 kg
Neck Diameter	120 mm	216 mm	216 mm	216 mm
Roller Base	O	O	O	O
Total 2ml Cryovial Capa.	750 ea	3000 ea	4800 ea	6000 ea

Roller Base

Cat. No.	Model	Description	EA Price ₩
C25-07-520	9902112	for BioRack 750	450,800
C25-07-530	9902147	for BioRack 3000, 4800, 6000	450,800

Other Accessory

Cat. No.	Model	Description	적용	EA Price ₩
C25-07-560	9702000	Cryoguard low level alarm	BioRack series	2,642,300
C25-07-570	9701321	Temperature logger		3,767,200

- Temperature logger alarm has an 3-way audible alarm for :
Low LN₂ Level, High Temperature, Low Battery Power



Shipping-type Cryosystem 랙 타입 액체 질소 저장 탱크

Cat. No.	Model	용량	Size	EA Price ₩
C25-07-250	Biotrek 3	4.2 L	Ø216×h455 mm	3,624,800
C25-07-260	Biotrek 10	13 L	Ø408×h658 mm	8,465,100

- Ultra Cold Condition 안에서 Sample의 이송용으로 적합함
- 내부에 Cryo-absorbant Material로 되어져 LN₂가 흐르지 않음
- 뚜껑을 양쪽에서 고정할 수 있음
- IATA Shipping Regulation 202, MDD Approved 93/42/EC망
- Statebourne / BioRack®

Cat. No.	C25-07-250	C25-07-260
Model	Biotrek 3	Biotrek 10
No. of Canister	1 ea	5 ea
용량	4.2 L	13 L
Evaporation Rate	0.15 %/day	0.7 %/day
Field Working Time	18 day	11 day
Neck Diameter	51 mm	216 mm
Weight Empty	5.1 kg	16.5 kg
Weight Full	8.5 kg	27 kg
Vial 용량 (Canister)	30 ea	510 ea
Vial 용량 (Cryobox)	-	500 ea

