

Micrometer Applicator 마이크로 어플리케이터

Cat. No.	Model	Film Width	EA Price ₩
A03-41-160	1117/100	100 mm	759,100
A03-41-168	1117/150	150 mm	759,100
A03-41-176	1117/200	200 mm	759,100
A03-41-184	1117/250	250 mm	843,400
A03-41-192	1117/300	300 mm	927,800


- 2 Micrometric Screws (From 0 ~ 13 mm, in 10 um Increment)
- 필름의 표면을 균일하게 코팅할때 사용함
- 요청에 의해 규격외 제품도 공급 가능
- LK Labkorea



Baker Applicator 베이커 어플리케이터

Cat. No.	Model	Film Width	EA Price ₩
A03-178-501	YBA-1	50 mm	939,800
A03-178-508	YBA-2	75 mm	939,800
A03-178-515	YBA-3	100 mm	939,800
A03-178-522	YBA-4	150 mm	1,096,400
A03-178-529	YBA-5	200 mm	1,253,100
A03-178-536	YBA-6	250 mm	1,409,700
A03-178-543	YBA-7	300 mm	1,566,300


- 도포의 두께를 조절 할 수 있는 가변형 Type의 Applicator
- 코팅두께 : 12.5 um, 25 um, 50 um, 75 um, 100 um, 125 um, 150 um, 175 um, 200 um, 225 um, 237.5 um, 250 um
- Yoshimitsu



4-Side Applicator, Sheen 4면 어플리케이터

Cat. No.	Model	Film Width	Gap Size	EA Price ₩
A03-41-112	1107/60/1	60 mm	30, 60, 90, 120 um	908,500
A03-41-117	1107/80/1	80 mm		939,800
A03-41-121	1107/60/2	60 mm	50, 100, 150, 200 um	908,500
A03-41-125	1107/80/2	80 mm		939,800

- Advantage of Combining Case of Use with 4 Gap Size in One unit
- By Simple Rotating The Applicator 90° The Next Gap is Placed on to The Test Surface
- Sheen



Vacuum Bed - Vacuum Plate 진공 베드

Cat. No.	Model	Bed Size (mm)	Coating Size (mm)	EA Price ₩
A03-41-081	VP-4230	w490×d300	w400×d250	1,018,100

- 재질 : Aluminum Anodizing
- Plate 위에 미세한 흡을 가지고 있어 Pump를 이용하여 압착
- LK Labkorea

