



Canister-type Cryosystem 캐니스터 타입 액체 질소 저장 탱크

- Strong, lightweight aluminium construction
- Statebourne

Extended Hold-time Model

Cat. No.	Model	용량	Size	EA Price ₩
C25-07-095	Bio 2	2 L	Ø185×h390 mm	2,293,100
C25-07-100	Bio 3	3 L	Ø216×h464 mm	2,536,300
C25-07-105	Bio 8	8 L	Ø407×h475 mm	3,000,600
C25-07-110	Bio 10	10 L	Ø305×h640 mm	3,124,100
C25-07-115	Bio 20	20 L	Ø407×h640 mm	3,512,400
C25-07-120	Bio 35	35 L	Ø480×h730 mm	4,350,500



- With narrow neck openings for high thermal efficiency
- These models offer long term storage possibilities with minimal LN₂ consumption

Cat. No.	C25-07-095	C25-07-100	C25-07-105	C25-07-110	C25-07-115	C25-07-120
Model	Bio 2	Bio 3	Bio 8	Bio 10	Bio 20	Bio 35
용량	2 L	3 L	8 L	10 L	20 L	35 L
Evaporation Rate	0.08 L/day	0.11 L/day	0.10 L/day	0.10 L/day	0.10 L/day	0.10 L/day
Static Hoding Time	25 day	27 day	80 day	100 day	200 day	350 day
Weight Empty	3 kg	4.3 kg	8.6 kg	7.9 kg	11.1 kg	18.6 kg
Weight Full	4.6 kg	6.7 kg	15.1 kg	16 kg	27.3 kg	46.9 kg
Neck Diameter	32 mm	51 mm	51 mm	51 mm	51 mm	51 mm
Roller Base	×	×	0	×	0	0
Number of Canister	3	6	6	6	6	6
Canister Dia. / Height	25 / 120 mm	38 / 120 mm	38 / 120 mm	38 / 275 mm	38 / 275 mm	38 / 275 mm
2ml Vial 용량	×	×	×	180 ea	180 ea	180 ea
Straw Capa. 0.25 / 0.5cc	360 / 120 ea	1200 / 500 ea	1200 / 500 ea	2400 / 1000 ea	2400 / 1000 ea	2400 / 1000 ea

Large 용량 Model

Cat. No.	Model	용량	Size	EA Price ₩
C25-07-125	Bio 12	12 L	Ø408×h470 mm	3,730,100
C25-07-130	Bio 21	21 L	Ø480×h579 mm	4,582,800
C25-07-135	Bio 22		Ø480×h537 mm	4,582,800
C25-07-140	Bio 34	34 L	Ø480×h710 mm	4,800,500
C25-07-145	Bio 36	35 L		5,497,000



- With wide neck openings for easy access
- These models allow the storage of large quantities of biological material

Cat. No.	C25-07-125	C25-07-130	C25-07-135	C25-07-140	C25-07-145
Model	Bio 12	Bio 21	Bio 22	Bio 34	Bio 36
용량	12 L	21 L	21 L	34 L	35 L
Evaporation Rate	0.19 L/day	0.19 L/day	0.28 L/day	0.18 L/day	0.30 L/day
Static Hoding Time	63 day	110 day	75 day	189 day	117 day
Weight Empty	10.8 kg	15.6 kg	15.7 kg	17.8 kg	19.1 kg
Weight Full	20.4 kg	32.6 kg	32.7 kg	45.3 kg	47.4 kg
Neck Diameter	70 mm	89 mm	120 mm	89 mm	120 mm
Roller Base	0	0	0	0	0
Number of Canister	9	6	10	6	10
Canister Dia. / Height	46 / 120mm	68 / 120mm	73 / 120mm	69 / 275mm	73 / 275mm
2ml Vial 용량	×	×	×	684 ea	1020 ea
Straw Capa. 0.25 / 0.5cc	4680 / 2200 ea	4920 / 2300 ea	8200 / 3800 ea	9840 / 4200 ea	16400 / 7800 ea

Roller Base

Cat. No.	Model	규격	EA Price ₩
C25-07-500	9902082	for Bio 8, 20, 12	375,700
C25-07-520	9902112	for Bio 21, 22, 34, 35, 36	450,800

Other Accessory

Cat. No.	Model	규격	EA Price ₩
C25-07-560	9702000	Cryoguard Low Level Alarm	2,642,300
C25-07-570	9701321	Temperature Logger	3,767,200



Roller Base