

## 1°C Cooler for 1.0 / 1.8ml Cryovial 1.0 / 1.8ml Cryovial용 / 1°C 쿨러

Cat. No.	Model	Hole	용량	Size (Ø×h)	EA Price ₩
C25-46-001	T525000	18	1.0 / 1.8 ml	117×86 mm	57,000

- 재질 : PC Body, HDPE Cover & Vial Holder
- Cell Cryo Preservation 및 Recovery용
- Cooling Rate : -1°C/min
- 100% IPA와 Mechanical Freeze가 필요함
- Tarsons



## Quick Cooler for Microcentrifuge Tube 퀵 쿨러, 0.5 ~ 1.5ml

Cat. No.	Model	Hole	용량	Size (w×d×h)	EA Price ₩
C25-46-003	T525040	12	0.2~1.5 ml	151×108×125 mm	94,000

- 재질 : PC Insulation Solution
- Specimer의 급속 냉각용 (DNA/RNA Sample)
- Dry Ice Solvent Bath의 대체용
- Cooling Temp : Down to -20°C, 최대 45분까지 유지 가능
- Tarsons



## -20°C Mini Cooler with Non Toxic Gel -20°C 미니 쿨러

- 재질 : PC, Non Toxic Gel
- Enzyme 및 Reagent들의 저온 유지용 및 Ice Bucket의 대체용
- Protect Critical Sample from Temperature Fluktuation in the Freezers
- 사용 방법 : -20 ~ -25°C 에서 24시간 이상 보관 후 사용
- Tarsons

### With Clear Cover, 기본형 -20°C 미니 쿨러

Cat. No.	Model	Hole	용량	Size (w×d×h, mm)	EA Price ₩
C25-46-016	T525030	12	1.5 ml	151×108×125	90,600
C25-46-018	T526030	32		243×157×146	165,700
C25-46-020	T526040	20	0.2 ml, PCR	151×108×125	146,100
C25-46-022	T525110	96			235,200

### With Gel Filled Cover, 겔 커버형 -20°C 쿨러

Cat. No.	Model	Hole	용량	Size (w×d×h, mm)	EA Price ₩
C25-46-024	T525060	32	1.5 ml	243×157×146	190,600



## 0°C Mini Cooler with Non Toxic Gel 0°C 미니 쿨러

- 재질 : PC, Non Toxic Gel
- Enzyme 및 Reagent들의 저온 유지용 및 Ice Bucket의 대체용
- Protect Critical Sample from Temperature Fluktuation in the Freezers
- 사용 방법 : -5 ~ -10°C 에서 24시간 이상 보관 후 사용
- Tarsons


### with Clear Cover, 기본형 0°C 미니 쿨러

Cat. No.	Model	Hole	용량	Size (w×d×h, mm)	EA Price ₩
C25-46-006	T525010	12	1.5 ml	151×108×125	90,600
C25-46-007	T526010	32		243×157×146	165,700
C25-46-008	T525020	12	15 ml	197×140×190	135,300
C25-46-009	T526050		Ø13 mm Tube		131,600

- 526050은 Blood Collection Tube와 같이 사용 가능

### with Gel Filled Cover, 겔 커버형 0°C 쿨러

Cat. No.	Model	Hole	용량	Size (w×d×h, mm)	EA Price ₩
C25-46-011	T525050	32	1.5 ml	243×157×146	180,600



- PCR Plate용 Cooler는 844 Page 참조
- Ice-free Cooling Box는 177 Page 참조

