




Measuring Cylinder, Amber Graduation BLAUBRAND®

메스 실린더, 갈색 눈금 Class A + Batch 보증서

Cat. No.	Model	용량	눈금	Tol.	높이	EA Price ₩	PK Price ₩
C27-59-101	32705	5 ml	0.1 ml	±0.05 ml	115 mm	14,900	2ea 27,700
C27-59-104	32708	10 ml	0.2 ml	±0.10 ml	140 mm	11,100	2ea 20,700
C27-59-107	32720	25 ml	0.5 ml	±0.25 ml	170 mm	12,000	2ea 22,400
C27-59-110	32728	50 ml	1.0 ml	±0.50 ml	200 mm	13,800	2ea 25,800
C27-59-113	32738	100 ml			260 mm	15,200	2ea 28,300
C27-59-116	32748	250 ml	2.0 ml	±1.00 ml	335 mm	23,700	2ea 44,300
C27-59-119	32754	500 ml	5.0 ml	±2.5 ml	365 mm	32,500	2ea 60,700
C27-59-122	32762	1 L	10.0 ml	±5.00 ml	465 mm	40,000	-
C27-59-125	32764	2 L	20.0 ml	±10.0 ml	505 mm	76,900	-

- Borosilicate Glass 3.3
- Complied with DIN EN ISO 4788, Calibrated to Contain (TC, In)
- 보증서 (Batch Certificate) 포함
- 갈색 (Amber) 침투 눈금 : 강한 세정 작업에도 갈색 눈금 유지
- With Spout and Hexagonal Base
- Brand GmBH




Measuring Cylinder, Blue Graduation BLAUBRAND®

메스 실린더, 파란 눈금 Class A + Batch 보증서

Cat. No.	Model	용량	눈금	Tol.	높이	EA Price ₩	PK Price ₩
C27-59-131	32105	5 ml	0.1 ml	±0.05 ml	115 mm	15,800	2ea 29,400
C27-59-134	32108	10 ml	0.2 ml	±0.10 ml	140 mm	11,900	2ea 22,200
C27-59-137	32120	25 ml	0.5 ml	±0.25 ml	170 mm	12,300	2ea 23,000
C27-59-140	32128	50 ml	1 ml	±0.5 ml	200 mm	14,300	2ea 26,700
C27-59-143	32138	100 ml			260 mm	15,600	2ea 29,000
C27-59-146	32148	250 ml	2 ml	±1.0 ml	335 mm	24,500	2ea 45,700
C27-59-149	32154	500 ml	5 ml	±2.5 ml	365 mm	33,400	2ea 62,400
C27-59-152	32162	1 L	10 ml	±5.0 ml	465 mm	40,000	-
C27-59-155	32164	2 L	20 ml	±10.0 ml	505 mm	76,900	-

- Borosilicate Glass 3.3
- Complied with DIN EN ISO 4788, Calibrated to Contain (TC, In)
- 보증서 (Batch Certificate) 포함
- 파란색(Blue) 침투눈금 : 강한 세정에도 파란 눈금 유지
- With Spout and Hexagonal Base
- Brand GmBH



Measuring Cylinder, USP Grade, BLAUBRAND®

메스 실린더, Class A + USP Batch 보증서

Cat. No.	Model	용량	눈금	Tol.	높이	EA Price ₩	PK Price ₩
C27-59-890	32805	5 ml	0.1 ml	±0.05 ml	115 mm	18,800	2ea 35,100
C27-59-893	32808	10 ml	0.2 ml	±0.10 ml	140 mm	18,800	2ea 35,100
C27-59-896	32820	25 ml	0.5 ml	±0.17 ml	170 mm	19,700	2ea 36,700
C27-59-899	32828	50 ml	1 ml	±0.25 ml	200 mm	21,200	2ea 39,600
C27-59-902	32838	100 ml		±0.5 ml	260 mm	23,000	2ea 42,900
C27-59-905	32848	250 ml	2 ml	±1.0 ml	335 mm	43,800	2ea 81,800
C27-59-908	32854	500 ml	5 ml	±2.0 ml	365 mm	74,600	2ea 139,300
C27-59-911	32862	1 L	10 ml	±3.0 ml	465 mm	108,700	-
C27-59-914	32864	2 L	20 ml	±6.0 ml	505 mm	184,300	-

- Borosilicate Glass 3.3
- Complied with DIN EN ISO 4788, Calibrated to Contain (TC, In)
- USP 보증서 (Batch Certificate) 포함
- With Spout and Hexagonal Base
- Brand GmBH

