


PP Low Form Measuring Cylinder PP 단형 실린더, Embossed Scale

Cat. No.	Model	용량	눈금	Tol.	높이	EA Price ₩	PK Price ₩
C27-59-807	41620	25 ml	0.5 ml	±0.5 ml	122 mm	4,400	10ea 40,500
C27-59-810	41628	50 ml	1.0 ml	±1.0 ml	142 mm	5,900	10ea 55,200
C27-59-813	41638	100 ml	2.0 ml	±2.0 ml	163 mm	6,600	10ea 61,600
C27-59-816	41648	250 ml	5.0 ml	±5.0 ml	192 mm	10,300	5ea 48,100
C27-59-819	41654	500 ml	10.0 ml	±10.0 ml	218 mm	13,800	5ea 64,200
C27-59-822	41662	1 L	20.0 ml	±20.0 ml	285 mm	19,200	5ea 89,500

- 재질 : PP (반투명), -10 ~ 135°C 내열
- Complied with DIN 12681 & ISO 6706, Calibrated to Contain (TC, In)
- 권장 사용 온도 : Max.80°C
- 권장 세정 온도 : Max.60°C (마크와 눈금 보호)
- Brand GmBH



PMP Measuring Cylinder PMP 메스 실린더, Class B



- 재질 : PMP (투명), 0 ~ 150°C 내열
- Complied with DIN 12681 & ISO 6706, Calibrated to Contain (TC, In)
- 121°C의 고압 멸균 후에도 정밀도 유지 • 권장 세정 온도 : Max.60°C (마크와 눈금 보호)
- Brand GmBH

Blue Scale, 청색 눈금

Cat. No.	Model	용량	눈금	Tol.	높이	EA Price ₩	PK Price ₩
C27-59-541	34708	10 ml	0.2 ml	±0.2 ml	145 mm	9,900	10ea 92,000
C27-59-543	34720	25 ml	0.5 ml	±0.5 ml	170 mm	12,800	10ea 119,400
C27-59-545	34728	50 ml	1.0 ml	±1.0 ml	200 mm	15,400	10ea 143,900
C27-59-547	34738	100 ml			250 mm	19,700	10ea 184,000
C27-59-549	34748	250 ml	2.0 ml	±2.0 ml	315 mm	29,500	5ea 137,600
C27-59-551	34754	500 ml	5.0 ml	±5.0 ml	360 mm	43,200	5ea 201,700
C27-59-553	34762	1 L	10.0 ml	±10.0 ml	440 mm	65,400	5ea 305,400
C27-59-555	34764	2 L	20.0 ml	±20.0 ml	535 mm	91,100	-

Embossed Scale, 몰드 눈금

Cat. No.	Model	용량	눈금	Tol.	높이	EA Price ₩	PK Price ₩
C27-59-561	34908	10 ml	0.2 ml	±0.2 ml	145 mm	6,800	10ea 63,100
C27-59-563	34920	25 ml	0.5 ml	±0.5 ml	170 mm	8,900	10ea 83,000
C27-59-565	34928	50 ml	1.0 ml	±1.0 ml	200 mm	10,700	10ea 100,000
C27-59-567	34938	100 ml			250 mm	13,700	10ea 127,700
C27-59-569	34948	250 ml	2.0 ml	±2.0 ml	315 mm	21,300	5ea 99,400
C27-59-571	34954	500 ml	5.0 ml	±5.0 ml	360 mm	32,400	5ea 151,300
C27-59-573	34962	1 L	10.0 ml	±10.0 ml	440 mm	49,800	5ea 232,400
C27-59-575	34964	2 L	20.0 ml	±20.0 ml	535 mm	77,100	-



PMP Measuring Cylinder PMP 메스 실린더, Class B

Cat. No.	Model	용량	눈금	Tol.	Size, $\Phi \times h$	EA Price ₩	PK Price ₩
C27-46-001	346010	10 ml	0.2 ml	±0.2 ml	13.6 x 147.2 mm	5,500	24ea 122,500
C27-46-003	346020	25 ml	0.5 ml	±0.5 ml	20.2 x 170.0 mm	5,300	18ea 91,100
C27-46-005	346030	50 ml	1.0 ml	±1.0 ml	25.8 x 199.5 mm	6,000	18ea 101,900
C27-46-007	346040	100 ml			30.9 x 249.5 mm	6,700	12ea 77,000
C27-46-009	346050	250 ml	2.0 ml	±2.0 ml	40.2 x 324.9 mm	11,200	12ea 128,000
C27-46-011	346060	500 ml	5.0 ml	±5.0 ml	52.2 x 388.0 mm	15,900	8ea 121,100
C27-46-013	346070	1 L	10.0 ml	±10.0 ml	65.5 x 458.1 mm	24,800	6ea 141,200
C27-46-015	346080	2 L	20.0 ml	±20.0 ml	83.4 x 533.7 mm	50,000	4ea 190,000

- 재질 : PMP (투명), 0 ~ 150°C 내열, Autoclavable
- Class B, Complied with DIN 12681 & ISO 6706, US FDA 21 CFR
- with Spout and Hexagonal Base
- Tarsons

