



PMP Measuring Cylinder, Blue Scale PMP 메스 실린더, Class A + Batch 보증서

Cat. No.	Model	용량	눈금	Tol.	높이	EA Price ₩	PK Price ₩
C27-59-581	35108	10 ml	0.2 ml	±0.10 ml	145 mm	15,300	2ea 28,600
C27-59-583	35120	25 ml	0.5 ml	±0.25 ml	170 mm	16,400	2ea 30,500
C27-59-585	35128	50 ml	1.0 ml	±0.5 ml	200 mm	17,200	2ea 32,000
C27-59-587	35138	100 ml			250 mm	23,400	2ea 43,700
C27-59-589	35148	250 ml	2.0 ml	±1.0 ml	315 mm	36,200	2ea 67,500
C27-59-591	35154	500 ml	5.0 ml	±2.5 ml	360 mm	54,800	2ea 102,400
C27-59-593	35162	1 L	10.0 ml	±5.0 ml	440 mm	73,000	-
C27-59-595	35164	2 L	20.0 ml	±10.0 ml	535 mm	120,200	-

- 재질 : PMP (투명), 0 ~ 150°C 내열
- Complied with DIN 12681 & ISO 6706, Calibrated to Contain (TC, In)
- 보증서 (Batch Certificate) 포함
- 121°C의 고압 멸균 후에도 정밀도 유지 • 권장 세정 온도 : Max.60°C (마크와 눈금 보호)
- Brand GmBH



PMP Measuring Cylinder PMP 메스 실린더, Class A

Cat. No.	Model	용량	눈금	Tol.	Size, Ø x h	EA Price ₩
C27-46-053	343020	25 ml	0.5 ml	±0.5 ml	20.2 x 170.0 mm	9,500
C27-46-056	343030	50 ml	1.0 ml	±1.0 ml	25.8 x 199.5 mm	10,100
C27-46-059	343040	100 ml			30.9 x 249.5 mm	13,900
C27-46-062	343050	250 ml	2.0 ml	±2.0 ml	40.2 x 324.9 mm	21,700
C27-46-065	343060	500 ml	5.0 ml	±5.0 ml	52.2 x 388.0 mm	32,900
C27-46-068	343070	1 L	10.0 ml	±10.0 ml	65.5 x 458.1 mm	47,000
C27-46-071	343080	2 L	20.0 ml	±20.0 ml	83.4 x 533.7 mm	100,200

- 재질 : PMP (투명), 0 ~ 150°C 내열, Autoclavable
- Class A, Complied with DIN 12681 & ISO 6706, US FDA 21 CFR
- with Spout and Hexagonal Base
- Tarsons



PFA Measuring Cylinder PFA 테프론 메스 실린더

Cat. No.	Model	용량	눈금	Size (Ø x h)	EA Price ₩
C27-64-221	F1119-01	5 ml	0.2 ml	8 x 140 mm	224,100
C27-64-223	F1119-02	10 ml		10 x 180 mm	259,600
C27-64-225	F1119-03	15 ml		12 x 190 mm	269,900
C27-64-227	F1119-04	20 ml		12 x 210 mm	349,600
C27-64-229	F1119-05	30 ml	0.5 ml	14 x 240 mm	393,700
C27-64-231	F1119-06	50 ml		19 x 240 mm	426,000
C27-64-233	F1119-07	100 ml	1.0 ml	26 x 250 mm	516,000
C27-64-235	F1119-08	200 ml	2.0 ml	33 x 290 mm	628,000
C27-64-237	F1119-09	300 ml	5.0 ml	41 x 290 mm	718,000
C27-64-239	F1119-10	500 ml		49 x 340 mm	852,000
C27-64-241	F1119-11	1 L		61 x 420 mm	1,412,200

- 재질 : PFA, -200 ~ 250°C 내열 • 최상의 내화학성을 지닌 반투명 테프론 재질
- Flon Industry



PTFE Measuring Cylinder PTFE 메스 실린더

Cat. No.	Model	용량	눈금	Size (Ø x h)	EA Price ₩
C27-460-001	A164-12	50 ml	1.0 ml	24 x 172 mm	638,700
C27-460-004	A164-16	100 ml	2.0 ml	30 x 202 mm	851,700
C27-460-007	A164-20	250 ml	2.5 ml	40 x 268 mm	1,184,900

- 재질 : PTFE & TFM, -200 °C to 250°C 내열
- with Spout, Reinforcing Ring and Hexagonal Base
- Sample은 단지 PTFE에만 노출되어짐
- Complied with ISO 4788
- BOLA

