

Glass Desiccator 유리 데시케이터

- 재질 : Sodalime Glass • Porcelain 중판 포함
- LK Labkorea

Clear Desiccator

Cat. No.	Model	상부 내경	전체 높이	중판 사이즈	EA Price ₩
D01-96-008	DC0150	Ø150 mm	h220 mm	Ø135 mm	20,700
D01-96-016	DC0180	Ø180 mm	h280 mm	Ø150 mm	26,700
D01-96-024	DC0210	Ø210 mm	h320 mm	Ø185 mm	39,200
D01-96-032	DC0240	Ø240 mm	h360 mm	Ø215 mm	60,000
D01-96-040	DC0300	Ø300 mm	h365 mm	Ø275 mm	90,600
D01-96-048	DC0360	Ø350 mm	h405 mm	Ø315 mm	132,900
D01-96-050	DC040	Ø400 mm	h470 mm	Ø355 mm	317,200

Amber Desiccator

Cat. No.	Model	상부 내경	전체 높이	중판 사이즈	EA Price ₩
D01-96-054	DA0150	Ø150 mm	h220 mm	Ø135 mm	36,100
D01-96-056	DA0180	Ø180 mm	h280 mm	Ø150 mm	44,900
D01-96-064	DA0210	Ø210 mm	h320 mm	Ø185 mm	64,200
D01-96-072	DA0240	Ø240 mm	h360 mm	Ø215 mm	101,400
D01-96-080	DA0300	Ø300 mm	h365 mm	Ø275 mm	154,500
D01-96-088	DA0360	Ø350 mm	h405 mm	Ø315 mm	189,800
D01-96-090	DA0360	Ø400 mm	h470 mm	Ø355 mm	518,800



Glass Vacuum Desiccator 유리 진공 데시케이터

- 재질 : Sodalime Glass
- Porcelain 중판, Glass Vacuum Valve 포함
- LK Labkorea

Clear Desiccator


Cat. No.	Model	상부 내경	전체 높이	중판 사이즈	EA Price ₩
D01-96-094	DC0150V	Ø150 mm	h220 mm	Ø135 mm	39,000
D01-96-104	DC0210V	Ø210 mm	h320 mm	Ø185 mm	66,000
D01-96-112	DC0240V	Ø240 mm	h360 mm	Ø215 mm	100,000
D01-96-120	DC0300V	Ø300 mm	h450 mm	Ø275 mm	150,000

Amber Desiccator

Cat. No.	Model	상부 내경	전체 높이	중판 사이즈	EA Price ₩
D01-96-134	DA0150V	Ø150 mm	h220 mm	Ø135 mm	57,800
D01-96-144	DA0210V	Ø210 mm	h320 mm	Ø185 mm	93,000
D01-96-152	DA0240V	Ø240 mm	h360 mm	Ø215 mm	147,600
D01-96-160	DA0300V	Ø300 mm	h450 mm	Ø275 mm	214,000

Optional PTFE Needle Valve

Cat. No.	Model	Description	EA Price ₩
D01-83-431	2822404	PTFE Needle Valve	54,800



Desiccator Plate 데시케이터 중판

Cat. No.	Model	Desiccator Size	Plate Size	EA Price ₩
D01-96-184	DP0180	Ø180 mm 용	150 mm	6,300
D01-96-192	DP0210	Ø210 mm 용	185 mm	7,000
D01-96-200	DP0240	Ø240 mm 용	215 mm	12,000
D01-96-208	DP0300	Ø300 mm 용	275 mm	15,700
D01-96-216	DP0360	Ø350 mm 용	315 mm	19,500

- 재질 : Porcelain
- LK Labkorea

