

## Desiccator, PP / PC 플라스틱 데시케이터

Cat. No.	Model	Plate Size	하부 id	전체 높이	EA Price ₩
D01-46-007	T401010	138 mm	150 mm	205 mm	50,100
D01-46-014	T401020	188 mm	198 mm	255 mm	94,800
D01-46-021	T401030	229 mm	240 mm	310 mm	131,900
D01-46-028	T401040	290 mm	296 mm	383 mm	199,900

- 재질 : PP 하부, Plate (White), PC Cover (투명)
- Large Interior Volume
- Thick Walled Desiccator
- Tarsons



## Vacuum Desiccator, PP/PC 플라스틱 진공 데시케이터

Cat. No.	Model	Plate Size	하부 id	전체 높이	EA Price ₩
D01-46-035	T402010	138 mm	150 mm	205 mm	83,100
D01-46-042	T402020	188 mm	198 mm	255 mm	110,800
D01-46-049	T402030	229 mm	240 mm	310 mm	131,900
D01-46-056	T402040	290 mm	296 mm	383 mm	256,600

- 재질 : PP 하부, Plate (White), PC Cover (투명)
- 3-way PP Stopcock, PTFE Key
- Red Silicon Gasket 포함
- Vacuum Connection : id6.4 mm Tubing용
- Hold a full Vacuum (740 mmHg) for 24 hr without Grease
- Tarsons



## Vacuum Desiccator, All Clear PC 투명 플라스틱 진공 데시케이터

Cat. No.	Model	Plate Size	하부 id	전체 높이	EA Price ₩
D01-46-063	T403010	138 mm	150 mm	205 mm	109,600
D01-46-070	T403020	188 mm	198 mm	255 mm	143,400
D01-46-077	T403030	229 mm	240 mm	310 mm	194,000
D01-46-084	T403040	290 mm	296 mm	383 mm	301,200

- 재질: PC (투명), PP Plate (White)
- 3-way PP Stopcock, PTFE Key
- Red Silicon Gasket 포함
- Vacuum Connection : id6.4 mm Tubing용
- Hold a Full Vacuum (740 mmHg) for 24 hr without Grease
- Tarsons



## Plastic Desiccator Accessory 플라스틱 데시케이터 악세서리

Cat. No.	Model	Description	EA Price ₩
D01-46-191	T401011	150 mm 용 od138 mm	10,000
D01-46-198	T401021	200 mm 용 od188 mm	18,800
D01-46-205	T401031	250 mm 용 od229 mm	25,900
D01-46-212	T401041	300 mm 용 od290 mm	39,500
D01-46-250	T402001	3-way PP Stopcock with PTFE Key ②	49,400
D01-46-307	T402011	150 mm 용	12,400
D01-46-314	T402021	200 mm 용	18,200
D01-46-321	T402031	250 mm 용	22,800
D01-46-328	T402041	300 mm 용	29,300

- Tarsons

