

Steel Auto Desiccator Cabinet

철제 자동 습도 조절 데시케이터, 20 ~ 50% RH

- Long Term Storage에 이상적인 경제형 철제 데시케이터
- Steel 재질의 선반은 높낮이 조절 가능함
- Humidity Range : 20 ~ 50%
- Analog 습도 조절
- 디지털 온습도계 장착
- Totech

1단형

배송비 별도

Cat. No.	D01-580-065	D01-580-070	D01-580-075
Model	CSD-151-20	CSD-252-20	CSD-502-20
용량	145 L	252 L	510 L
Humidity Control	20 ~ 50% RH, Analog Control		
Size (w×d×h, mm)	내부	490×530×560	870×330×830
	외부	500×580×630	880×380×900
Size (w×d×h, mm)	870×690×830		
외부	880×740×900		
선반 수	3개		
Door 수	1개	2개	
온습도계	Digital 기본 장착		
Caster	X	X	X
재질	Steel Body, Glass Window		
Power	AC 220V		
인증	CE		
EA Price ₩	1,321,800	1,414,100	2,766,700



CSD-1106-20

CE

다단형

배송비 별도

Cat. No.	D01-580-080	D01-580-085
Model	CSD-302-20	CSD-702-20
용량	340 L	680 L
Humidity Control	20 ~ 50% RH, Analog Control	
Size (w×d×h, mm)	내부	490×590×1160
	외부	500×640×1230
Size (w×d×h, mm)	610×730×1560	
	620×780×1840	
선반 수	7개	5개
Door 수	2개	
온습도계	Digital 기본 장착	
Caster	X	O
재질	Steel Body, Glass Window	
Power	AC 220V	
인증	CE	
EA Price ₩	1,998,200	3,227,800



CSD-1104-20

CE



CSD-702-20

CE



CSD-302-20

CE

다단형

배송비 별도

Cat. No.	D01-580-090	D01-580-095
Model	CSD-1104-20	CSD-1106-20
용량	1160 L	1160 L
Humidity Control	20 ~ 50% RH, Analog Control	
Size (w×d×h, mm)	내부	1190×620×1560
	외부	1200×670×1840
선반 수	5개	
Door 수	4개	6개
온습도계	Digital 기본 장착	
Caster	O	O
재질	Steel Body, Glass Window	
Power	AC 220V	
인증	CE	
EA Price ₩	4,303,700	4,457,500



CSD-502-20

CE



CSD-252-20

CE



CSD-151-20

CE