

PFA Evaporating Dish PFA 테프론 증발 접시

Cat. No.	Model	용량	Size	EA Price ₩
D02-59-151	45800	25 ml	∅50×h25 mm	78,100
D02-59-155	45802	50 ml	∅50×h54 mm	103,700

- 재질 : PFA, -196 ~ 250°C 내열, Autoclavable
- Snap-On Lid 재질 : PE, Non-Autoclavable
- Brand GmBH



PTFE Evaporating Dish PTFE 테프론 증발 접시, Conical shape

Cat. No.	Model	용량	O.D. Top Size	O.D. Bottom Size	EA Price ₩
D02-460-031	A169-01	25 ml	∅ 38 x h 37mm	∅ 34mm	53,600
D02-460-035	A169-02	50 ml	∅ 49 x h 41mm	∅ 46mm	54,600
D02-460-039	A169-03	100 ml	∅ 64 x h 41mm	∅ 60mm	74,600
D02-460-043	A169-04	250 ml	∅ 97 x h 62mm	∅ 90mm	101,600

- 재질 : PTFE, -200 °C to 250°C 내열
- Conical Shape, with Spout
- BOLA



PTFE Evaporating Dish PTFE 테프론 증발 접시, Cylindrical shape

Cat. No.	Model	용량	O.D. Size	I.D. Size	EA Price ₩
D02-460-061	A176-02	100 ml	∅56×h18 mm	∅100 mm	87,000
D02-460-063	A176-03	250 ml	∅130×h36 mm	∅125 mm	169,600


- 재질 : PTFE, -200 °C to 250°C 내열
- Cylindrical Shape, with Spout
- BOLA



PTFE Evaporating Dish PTFE 테프론 증발 접시

Cat. No.	Model	용량	Size	EA Price ₩
D02-59-222	009.005	5 ml	∅27×h16 mm	14,200
D02-59-226	009.025	25 ml	∅44×h26 mm	22,000
D02-59-230	009.075	75 ml	∅56×h41 mm	28,600
D02-59-234	009.100	100 ml	∅60×h45 mm	40,700

- 재질 : PTFE, -200 ~ 280°C 내열
- 최상의 내화학적, 내열성, 내마모성
- Crucible Form Without Spout
- Cowie



Casserole with Handle, Porcelain 자제 카세롤

Cat. No.	Model	용량	Size	EA Price ₩
D02-128-411	209/1	22 ml	∅45×h26 mm	11,100
D02-128-415	209/2	84 ml	∅67×h39 mm	14,600
D02-128-419	209/3	176 ml	∅85×h49 mm	18,700
D02-128-423	209/4	270 ml	∅101×h51 mm	21,100
D02-128-427	209/5	450 ml	∅120×h64 mm	26,700
D02-128-431	209/6	960 ml	∅150×h86 mm	44,700

- Max.Temp. : 1,000°C
- 손잡이가 부착된 자제 카세롤
- Jipo / LabCeramic®

