

Stainless Steel Evaporating Dish 스테인레스 증발 접시

- 재질 : 18/10 Stainless Steel, Antimagnetic, Chemical-Resistant
- Usbeck

High Shape

Cat. No.	Model	용량	Size (Ø×h)	EA Price ₩
D02-284-117	5030	100 ml	80×40 mm	34,500
D02-284-121	5031	250 ml	100×47 mm	38,400
D02-284-125	5032	500 ml	120×57 mm	58,800
D02-284-129	5033	1 L	160×65 mm	93,400
D02-284-133	5034	2 L	200×90 mm	150,000

Glass Oil Bath (Dish) 유리 오일 배스

Cat. No.	Model	Capacity	Size (Ø×h)	EA Price ₩
B04-91-245	GB020	2 L	160×135 mm	120,000
B04-91-255	GB045	4.5 L	200×165 mm	170,000

- Borosilicate Glass 3.3, USP Type I, ASTM E438, Type I, Class A
- Heavy Wall Type의 고온용 유리 수조
- Digital Oil Bath의 부품으로 사용되거나 독립적으로 사용 가능 (80 Page 참조)
- LK Labkorea

Stainless Steel Dish 스테인레스 디쉬

Cat. No.	Model	Size	EA Price ₩
D02-132-101	1100-D-01	Ø100×h65 mm	170,800
D02-132-104	1100-D-02	Ø150×h75 mm	345,400
D02-132-107	1100-D-03	Ø200×h100 mm	498,800

- Stainless Steel
- Hotplate Stirrer에서 Oil Bath용으로 사용 가능
- Cover 포함
- Chemglass

Compartment Glass Petri Dish, Simax® 칸막이 유리 페트리 디쉬

Cat. No.	Model	Section	Size	EA Price ₩	PK Price ₩
D02-52-561	491006120	3 칸	Ø120×h20 mm	37,100	12ea 444,900
D02-52-569	491007080	4 칸	Ø80×h15 mm	31,700	24ea 710,200
D02-52-571	491007120		Ø120×h20 mm	39,700	12ea 444,900

- Simax® Borosilicate Glass 3.3, USP Type I, ASTM E438, Type I, Class A
- Heavy Wall Type으로 두껍고 충격에 강함
- 고품질의 내열, 내화학성 유리 (450°C 내열)
- KAVALIERRGLASS

Compartmented PS Petri Dish 칸막이 PS 페트리 디쉬

Cat. No.	Model	Section	Size	CS Price ₩
D02-491-071	752011	2칸×20 ml	Ø90×h15 mm	500ea 82,900
D02-491-081	752021	3칸×13 ml		500ea 99,300
D02-491-091	752031	4칸×10 ml		500ea 99,300

- Polystyrene 재질로 투명도가 높아 내용물 관찰이 쉬움
- 균일한 바닥으로 수평도가 우수함
- 전자빔 처리로 완벽한 살균 효과 (Non-Pyrogenic, DNase/RNase-free)
- LK Labkorea



Flat Form

Cat. No.	Model	용량	Size (Ø×h)	EA Price ₩
D02-284-101	5010	60 ml	60×30 mm	30,500
D02-284-105	5011	100 ml	75×30 mm	34,500
D02-284-109	5012	175 ml	100×30 mm	40,600
D02-284-113	5013	400 ml	150×26 mm	44,600

