

PFA Petri Dish PFA 테프론 페트리 디쉬

Cat. No.	Model	Size	EA Price ₩
D02-64-137	F1078-01	∅50×h14 mm	47,500
D02-64-140	F1078-02	∅75×h17 mm	83,100
D02-64-143	F1078-03	∅100×h20 mm	107,000


- 재질 : PFA, -200 ~ 260 °C 내열
- Flon Industry



PTFE Petri Dish with Cover PTFE 테프론 페트리 디쉬와 커버

Cat. No.	Model	Size	EA Price ₩
D02-64-237	F1080-01	∅100×h25 mm ①	385,200
D02-64-240	F1081-01	∅90×h60 mm ②	706,000

- 재질 : PTFE, -200 ~ 280°C 내열
- Flon Industry



Quartz Petri Dish 석영 페트리 디쉬

Cat. No.	Model	Size	EA Price ₩
D02-60-101	QPD06015	∅60×h15 mm	143,000
D02-60-104	QPD09017	∅90×h17 mm	195,000
D02-60-107	QPD12025	∅120×h25 mm	221,000
D02-60-111	QPD15030	∅150×h30 mm	260,000

- 재질 : Quartz
- 사용 온도 : Max. 1250°C
- LK Labkorea



Stainless Steel Petri Dish 스테인레스 페트리 디쉬 [샤레]

Cat. No.	Model	Size	EA Price ₩
D02-102-117	CYS1223	∅60×h20 mm	7,200
D02-102-126	CYS1224	∅75×h20 mm	8,400
D02-102-135	CYS1226	∅90×h20 mm	10,200
D02-102-144	CYS1227	∅120×h20 mm	15,700
D02-102-153	CYS1228	∅170×h20 mm	18,100

- 재질 : Stainless Steel 304
- LK Labkorea



Glass Conway Dish 유리 콘웨이 디쉬

Cat. No.	Model	Size (∅×h)	Chamber (∅×h)	EA Price ₩
D02-11-107	GCD60	60×8 mm	35×4 mm	44,500

- 재질 : Sodalime Glass
- 시료의 확산과 반응을 확인하는 디쉬
- LK Labkorea



Conway Diffusion Cell 콘웨이 디쉬

Cat. No.	Model	Size (∅×h)	Chamber (∅×h)	PK Price ₩
D02-11-101	40941-0000	83×4 mm	60×10 mm	3ea 72,700

- 재질 : PP
- 혈액이나 혈청안의 이산화탄소, 암모니아, 알콜 확산 측정용
- Bel-Art

