

"Dispensette® S TA" Bottle Top Dispenser 불산 및 고순도용 바들 탑 디스펜서

- Trace Analysis를 위한 High-Purity Acids, Alcalis or Hydrogen Peroxide 용
- Hydrofluoric acid (HF) 용으로 사용 (Platinum-Iridium Valve Spring)
- Analog Adjustment
- 내부에 Sealing Ring이 없어 유지 보수 및 세척이 간단함
- ISO9001과 GLP에 따른 보정 작업이 가능
- 2가지 Valve Spring
 - Platinum-Iridium : Not Suitable Hydrogen Peroxide
 - Tantal : Not Suitable Hydrofluoric Acid and Sodium hydroxide.30%
- 공통 구성품 : Dispenser with GL45 Thread, Conformity Certified, Performance Certificate, Discharging Tube, Telescoping Filling Tube, Mounting Tool, Adapter (GL25, GL28, GL32, GL38, S40 용)
- Not Autoclavable, Not Suitable for Organic Solvent
- Brand GmBH



Without Recirculation Valve

Cat. No.	Model	용량	Value Spring	눈금	Accuracy	Coefficient of Variation	EA Price ₩
D03-59-600	4640040	1 ~ 10 ml	Platin-Iridium	0.2 ml	≤ ±0.5% (50 ul)	≤ ±0.1% (10 ul)	1,469,200
D03-59-604	4640240	1 ~ 10 ml	Tantal	0.2 ml	≤ ±0.5% (50 ul)	≤ ±0.1% (10 ul)	1,469,200

With Recirculation Valve

Cat. No.	Model	용량	Value Spring	눈금	Accuracy	Coefficient of Variation	EA Price ₩
D03-59-650	4640041	1 ~ 10 ml	Platin-Iridium	0.2 ml	≤ ±0.5% (50 ul)	≤ ±0.1% (10 ul)	1,565,200
D03-59-654	4640241	1 ~ 10 ml	Tantal	0.2 ml	≤ ±0.5% (50 ul)	≤ ±0.1% (10 ul)	1,565,200

Dispensing Cartridge with Safety Ring ①

Cat. No.	Model	용량	보증서	EA Price ₩
D03-59-660	708035	10 ml	포함	874,900



Discharge Tube

with Recirculation Valve ②

Cat. No.	Model	Valve Spring	Length	EA Price ₩
D03-59-664	708122	Platinum-iridium	105 mm	189,400
D03-59-668	708124	Tantalum		189,400

without Recirculation Valve ③

Cat. No.	Model	Valve Spring	Length	EA Price ₩
D03-59-672	708022	Platinum-iridium	105 mm	31,100
D03-59-676	708024	Tantalum		31,100


Flexible Discharge Tube with Recirculation Valve ④

Cat. No.	Model	용량	Tube Size	EA Price ₩
D03-59-949	708132	10 ml	od3 × id2 mm	198,300

- PTFE Coiled, Length Approx. 800 mm, with Safety Handle

PTFE Threaded Adapter for GL45

Cat. No.	Model	Bottle Thread	Top Thread	EA Price ₩
D03-460-200	H978-30 ⑤	GL32	GL45	96,000
D03-460-205	H978-40 ⑥	GL40/S40		38,700



- PFA Bottle은 147 ~ 150 Page 참조

