

Amber Square Dispenser Bottle 갈색 사각 디스펜서용 용매병

- 재질 : Amber Soda-Lime Glass
- 2.5 L (704010)은 Round Shape (원형)임
- PP Screw Cap과 LDPE Pouring Ring 포함
- Brand GmBH

Coated

Cat. No.	Model	용량	Size	Thread	EA Price ₩
D03-59-901	704002	100 ml	∅50×h125 mm	GL32	21,600
D03-59-903	704004	250 ml	∅65×h160 mm		28,700
D03-59-905	704006	500 ml	∅80×h195 mm		29,900
D03-59-907	704008	1 L	∅95×h230 mm	GL 45	36,100
D03-59-909	704010	2.5 L	∅140×h300 mm		66,300

- Max. Operating Temp. (Coated) : 80°C
- Ethylene-Acrylate 코팅 제품은 제품이 파손되었을 때 사용자의 위험을 줄여줌
- 코팅 제품은 60°C 이상에서 세척시 코팅이 벗겨질 수 있음

Uncoated

Cat. No.	Model	용량	Size	Thread	EA Price ₩
D03-59-911	704012	100 ml	∅50×h125 mm	GL 32	11,400
D03-59-913	704014	250 ml	∅65×h160 mm		13,900
D03-59-915	704016	500 ml	∅80×h195 mm		14,500
D03-59-917	704018	1 L	∅95×h230 mm	GL 45	17,200
D03-59-919	704020	2.5 L	∅140×h300 mm		37,200



Thread Bottle Adapter for Dispenser 디스펜서용 바틀 어댑터

- 재질 : PP, ETFE, PTFE
- Outer Thread GL 32 or GL 45 (Dispenser 연결 규격)
- Brand GmBH

GL 32

Cat. No.	Model	재질	Bottle Thread/ Ground Joint	EA Price ₩
D03-59-865	704322	PP	GL 22	15,700
D03-59-867	704325		GL 24, 25	14,500
D03-59-869	704328		GL 28	14,500
D03-59-871	704330		GL 30	18,900
D03-59-873	704338		GL 38	21,100
D03-59-875	704340		GL 40	22,300
D03-59-879	704345		GL 45	22,900
D03-59-885	704375	ETFE	GL 24, 25	32,900
D03-59-887	704378		GL 28	32,900
D03-59-889	704388		GL 38	40,000
D03-59-891	704390		GL 40	42,900
D03-59-962	704395		GL 45	51,000
D03-59-976	704419	PP	NS 19/26	24,900
D03-59-978	704424		NS 24/29	26,900
D03-59-980	704429		NS 29/32	27,700

GL 45

Cat. No.	Model	재질	Bottle Thread	EA Price ₩
D03-59-877	704343	PP	S 40	22,300
D03-59-881	704348		S 60	64,600
D03-59-883	704350		S 50	55,200
D03-59-960	704391	PTFE	S 40	50,500
D03-59-964	704396	PP	GL 32-33	23,300
D03-59-966	704397		GL 38	23,300
D03-59-968	704398	ETFE	GL 32-33	45,200
D03-59-970	704399		GL 38	45,500
D03-59-972	704430	PP	S 54	54,700
D03-59-974	704431		GL 35	24,000

