

Universal HF Dispenser, POLYFIX®

불산 검용 바틀 탑 디스펜서

Cat. No.	Model	용량	눈금	EA Price ₩
D03-66-405	10110633	1 ~ 5 ml	0.1 ml	1,662,400
D03-66-411	10110638	2 ~ 10 ml	0.2 ml	1,662,400
D03-66-430	10110644	5 ~ 30 ml	0.5 ml	2,468,000
D03-66-435	10110648	10 ~ 50 ml	1.0 ml	3,701,800

- 2 Step의 용량 조절 Valve가 있어 분주 용량을 정밀하게 Setting
- 불산(Hydrofluoric Acid)을 포함한 모든 Liquid Chemical용 디스펜서
- PTFE Cylinder & Piston, Valve Block
- Accuracy : $\leq \pm 0.6\%$ (Tolerance : $\leq \pm 0.2\%$)
- Solid PTFE로 만들어진 GL32 Fix Adapter
- GL40 & GL45 PTFE Adapter
- Poulten & Graf



Optional Accessory for OPTIFIX®

OPTIFIX® 용 약세사리

- Poulten & Graf

Flexible Universal Drain Tube

Cat. No.	Model	적용	EA Price ₩
D03-66-923	10110401	2 ~ 10 ml용 Dispenser	241,800
D03-66-927	10110402	30 ~ 50 ml용 Dispenser	302,200
D03-66-931	10110406	100ml용 Dispenser	302,200

- 재질 : PTFE

Discharge Tube

Cat. No.	Model	적용	길이	EA Price ₩
D03-66-953	10109427	2 ml용 Dispenser	90 mm	62,400
D03-66-957	10109438	5 ~ 10 ml용 Dispenser		62,400
D03-66-961	10109448	30 ~ 50 ml용 Dispenser	99 mm	68,100
D03-66-964	10109451	100 ml용 Dispenser	116 mm	94,800
D03-66-967	10109458	200 ~ 500 ml용 Dispenser		135,300


Ground Joint Adapter, PP

Cat. No.	Model	적용	EA Price ₩
D03-66-982	10109224	ST 24	36,300
D03-66-986	10109229	ST 29	40,500
D03-66-990	10109245	ST 45	68,400

- GL32 와 Ground Joint의 연결용 어댑터



Flexible Universal Drain Tube



Discharge Tube



Adapter

Dispenser Bottle

디스펜서용 용매병, LukeGL®

Cat. No.	Model	용량	Center	Side	EA Price ₩
D03-11-555	DB0500	500 ml	GL45	GL14	31,600
D03-11-562	DB1000	1 L			42,800
D03-11-569	DB2000	2 L			56,300
D03-11-576	DB5000	5 L			196,900

- Borosilicate Glass 3.3
- Side Neck Size : GL14
- LukeGL®

