

## Handy Aspirator - Hand Vacuum Pump & Aspirator 핸디 진공 펌프 & 핸디 아스피레이터

- 가볍고 Field 용으로 편리한 Vacuum & Pressure Pump
- Vacuum Release Valve 포함
- Connection: id6 mm Tubing에 알맞음
- Mityvac®

### Vacuum Pump & Pressure Pump 겸용, 16 ml/stroke

Cat. No.	Model	Pumping Rate	Max. Vacuum	EA Price ₩
A04-53-030	MV-8255	16 ml/stroke	25 inHg (635 mmHg)	189,400

- Max. Pressure : 12 ~ 15 psi (Top Port)
- Vacuum & Pressure Gauge 기본 장착

### Vacuum Gauge only, 32 ml/stroke

Cat. No.	Model	Pumping Rate	Max. Vacuum	EA Price ₩
A04-53-052	MVA-6178	Vacuum Gauge, 0 ~100 kpa		72,800



## Plastic Aspirator, Water Jet pump 플라스틱 아스피레이터

Cat. No.	Model	재질	길이	EA Price ₩
A04-59-911	159600	Plastic	210 mm	78,600

- 재질 : PP, FKM, PTFE (High Chemical Resistance)
- Operating Temperature : Max.80 °C (Water)
- 역류 방지 밸브 (Non-return Valve) 장착
- Very Low Water Consumption : 190 L/h (3.5 bar 물 공급시)
- Constant Ultimate Pressure : 16 mbar (3 ~ 6 bar, 12°C 물 공급시)
- High Suction Capacity : Approx. 400 L/h of Air (3.5 bar, 12°C 물 공급시)
- 구성 어댑터 : Sleeve Nut R3/4", Reducing Adapter R1/2",  
od10 ~ 12 mm Hose Connector,  
od6 ~ 9 mm 분리형 Vacuum Connector (GL14)

• Brand GmBH

### Optional Reducing Adapter

Cat. No.	Model	Description	EA Price ₩
A04-59-913	159665	R3/4" to R3/8", PP	8,600
A04-59-915	159670	R3/4" to M22 Thread, PP	8,600



## Glass Aspirator 유리 아스피레이터

Cat. No.	Model	재질	Max. Vacuum	EA Price ₩
A04-83-222	5580003	Glass	Up to 16 mbar	183,600

- Borosilicate Glass 3.3, USP Type I, ASTM E438, Type I, Class A
- 16 mbar의 진공도까지 가능함
- Lenz

