

## Hample Distilling Column 햄플형 증류 컬럼, LukeGL®


- Borosilicate Glass 3.3
- 길이 : 상, 하 Joint부의 길이를 제외한 컬럼 길이
- LukeGL®

### 일반형 ①

Cat. No.	Model	길이	Joint	EA Price ₩
C17-83-008	5115200	100 mm	14/23	25,300
C17-83-016	5115201	150 mm	24/40	30,400
C17-83-024	5115202	200 mm		33,700
C17-83-032	5115203	300 mm		33,700

### 자켓형, Vacuum Jacketed ②


Cat. No.	Model	길이	Joint	EA Price ₩
C17-83-040	5119200	100 mm	14/23	42,200
C17-83-048	5119201	150 mm	24/40	47,200
C17-83-056	5119203	200 mm		50,600
C17-83-064	5119204	300 mm		53,700



## Distilling Column, Vacuum Jacketed 자켓형 증류 컬럼, LukeGL®

Cat. No.	Model	길이	Joint	EA Price ₩
C17-83-072	5115401	150 mm	14/23	59,000
C17-83-080	5115402	200 mm	24/40	62,400
C17-83-088	5115403	250 mm		75,900
C17-83-096	5115404	300 mm		87,700

- Borosilicate Glass 3.3
- Simple Design, Particularly Useful in Fractionations
- 길이 : 상, 하 Joint부의 길이를 제외한 컬럼 길이
- LukeGL®



## Kugelrohr Distilling Bulb Kugelrohr 증류 벌브, LukeGL®

Cat. No.	Model	용량	Joint	EA Price ₩
C17-83-104	L122710	10 ml	14/23	72,900
C17-83-112	L122712	25 ml		75,400
C17-83-120	L122713	50 ml		85,200
C17-83-128	L122714	100 ml		85,200

- LukeGL®

### Optional Accessory

Cat. No.	Model	Description	EA Price ₩
A01-83-088	5092001	14/23 Straight Inlet Adapter	11,100
F04-10-010	L110101E	10 ml, 14/23 1-N Round Flask	9,600
F04-10-012	L110102E	25 ml, 14/23 1-N Round Flask	10,800

## Snyder Distilling Columns 스나이더형 증류 컬럼, LukeGL®

Cat. No.	Model	용량	Joint	EA Price ₩
C17-83-136	5136701	10 ml	14/23	47,300
C17-83-144	5136702	25 ml	24/40	50,700
C17-83-152	5136703	50 ml		54,100
C17-83-160	5136704	100 ml		58,400

- 주로 Kuderna Danish Flask와 같이 사용

- LukeGL®

