

Distillation Receiver, Dean-Stark Trap 딘스타크 트랩, Moisture Test

Cat. No.	Model	용량	Joint	EA Price ₩
C17-132-768	1265-06	20 ml	24/40	745,200

- Borosilicate Glass 3.3, USP Type I, ASTM E438, Type I, Class A
- For Solvents that are either lighter a heavier than water
- Chemglass



Distillation Receiver, Dean-Stark Trap 딘스타크 트랩

Cat. No.	Model	Joint	Lower Reservoir	EA Price ₩
C17-132-786	1266-10	24/40	50 ml	822,200

- Borosilicate Glass 3.3, USP Type I, ASTM E438, Type I, Class A
- For azeotropic distillation using solvent, heavier than water
- 용량 (Stopcock) 밑의 용량임
- Chemglass



Distillation Receiver, Solvent Recovery 솔벤트 증류 리시버

Cat. No.	Model	용량	Joint	EA Price ₩
C17-132-792	1272-01	250 ml	24/40	1,045,300
C17-132-798	1272-02	500 ml		1,003,300

- Borosilicate Glass 3.3, USP Type I, ASTM E438, Type I, Class A
- Solvent의 제거 및 환원용
- Chemglass



Semi-Micro Kjeldal Quantitative Apparatus 질소 정량 장치 / 세미 마이크로 킬달장치, LukeGL®

Cat. No.	Model	Description	EA Price ₩
C17-132-828	1096-81	질소 정량 장치	204,800

- 대한 약전 규격
- Borosilicate Glass 3.3
- LukeGL®

Component Part

Cat. No.	Model	Part	Description	EA Price ₩
C17-132-834	1096-91	A	킬달 플라스크, 250 ml	31,300
C17-132-840	1096-92	B	동근 플라스크, 1 L	30,100
C17-132-846	1096-93	C	연결관	12,000
C17-132-852	1096-94	D	깔대기와 코크	24,100
C17-132-858	1096-95	E	냉각기, 300 mm	40,100
C17-132-864	1096-96	F	깔대기와 코크, 알칼리용액 투입용	60,700
C17-132-870	1096-97	G	삼각 플라스크, 250 ml	5,100

