


Sulfide Determination Apparatus 아황산 측정장치 / 이산화황 측정장치, LukeGL®

Cat. No.	Model	Description	EA Price ₩
C17-132-876	1096-61	1 L Complete Set	191,700

- 대한 약전 규격
- Borosilicate Glass 3.3
- LukeGL®

Component Part

Cat. No.	Model	Part	Description	EA Price ₩
C17-132-882	1096-62	A	Inlet Adapter 24/40	15,500
C17-132-888	1096-63	B	Addition Funnel, 250 ml	56,600
C17-132-894	1096-64	C	3-N Flask, 1 L, 24/40	59,000
C17-132-900	1096-65	D	Gas Inlet Adapter, 24/40	18,400
C17-132-906	1096-66	E	Allihn Condenser 300 mm, 24/40	39,900
C17-132-912	1096-67	F	Bubbler Tube, 24/40	20,200
C09-83-030	KC-24	G	Joint Clamp, 24/40	6,000
C27-59-113	32738	H	Cylinder, 100 ml	15,200



Fluorine Distilling Apparatus 불소 측정장치 / 수증기 측정장치, LukeGL®


- Borosilicate Glass 3.3
- LukeGL®

Cat. No.	Model	Description	EA Price ₩
C17-132-936	1096-71	1 L Complete Set	314,000

- 1L Complete Set에는 TH1420 2개 포함됨

Component Part

Cat. No.	Model	Part	Description	EA Price ₩
C17-132-952	1096-100	A	1-Neck Round Flask 1 L	27,800
C17-132-956	1096-101	B	상 아답타	60,000
C17-132-960	1096-102	C	킬달 아답타	39,400
C17-132-964	1096-103	D	냉각기	48,200
C17-132-968	1096-104	E	크라이젠 아답타	39,400
C17-132-972	1096-105	F	역류방지관	24,100
C17-132-976	1096-106	G	온도계 꽃이	24,100
A01-143-100	TH1420	G	테프론 온도계 어댑터	8,000
C17-132-980	1096-107	H	2-Neck Round Flask, 500 ml	28,700
C27-59-113	32738	I	Cylinder, 100 ml	15,200



- Distillation Cart System은 870 Page 참조

