

Multi-Function DO-BOD Meter 에서 계속...

General Accessory

Cat. No.	Model	Description	EA Price ₩
P08-162-011	30019823	uPlace™: Electrode Holder (Incl. in standard delivery)	
C19-162-904	30041154	Semi-transparent cover (Incl. in standard delivery)	
C19-162-916	30041155	Protective Film for Display (2 pcs)	
C19-162-902	11120270	Main Adapter 12V (Incl. in standard delivery)	
P08-162-020	30040000	uMix™ Magnetic Stirrer	
D06-162-521	51192107	LogStraight™ Fingerprint Reader	
C19-162-922	21901309	USB Cable for Barcode Reader	
C19-162-914	51302888	Rondolino Sample Changer	

Portable Dissolved Oxygen Meter 휴대용 DO 미터, Seven2Go S4/S9

Cat. No.	Model	형식	전극	EA Price ₩
D08-162-420	S4-Standard Kit	Polarographic, General	InLab® 605-ISM	
D08-162-430	S4-Field Kit	Polarographic, 현장 측정용		
D08-162-440	S9-Standard Kit	Optical, General	InLab® OptiOx™	
D08-162-450	S9-Field Kit	Optical, 현장 측정용		

- 측정 범위
 - DO : 0.000 ~ 99.9 mg/L (S4) / 0.000 ~ 50.0 mg/L (S9)
 - 온도 : 0.0 ~ 60.0°C (S4), 0.0 ~ 50.0°C (S9)
- 2가지 측정 방법 :
 - 폴라로그래픽 DO : 매우 복잡한 수용액 샘플에서도 DO 측정 가능 (S4)
 - 광학 DO : 매우 정확하고 신속한 InLab® OptiOx 센서를 지원 (S9)
- 자동 센서 인식 : 오류 없는 교정 및 센서 취급을 위해 모든 InLab® ISM 센서 지원
- IP67 등급 : 견고한 방수 하우징
- USB Interface를 이용한 데이터 전송 (S9)
- 한국어 지원 (S9)
- Mettler Toledo

Model		S4	S9
DO	Mesuring Range	0.00 ~ 99.99 mg/L (ppm)	0.00 ~ 50.00 mg/L (ppm)
	Resolution	0.01 mg/L	
	Accuracy with Sensor	±0.2 mg/L	±0.1 mg/L (0 ~ 8 mg/L) ±0.2 mg/L (8 ~ 20 mg/L) ±10% (20 ~ 50 mg/L)
Temperature	Measuring Range	0.0 ~ 60.0°C	0.0 ~ 50.0°C
	Resolution	0.1°C	0.1°C
	Accuracy	±0.2°C	±0.1°C
Calibration		2 Point	
Language		English	12 Language (incl. Korean)
측정 data 저장		200	2000 (GLP Conform)
Protecting Rating		Watertight and dustproof (IP67)	
Size		w222 × d70 × h35 mm	
Weight		270 g	290 g
Sensor Input		BNC / RCA (cinch)	Mini LTW (IP67)
Security and Compliance		-	Password protection GLP User management
Interface		No	Micro USB (USB Device)

Accessory

Cat. No.	Model	Description	EA Price ₩
C19-162-955	30122300	Seven2Go uGo 휴대용 케이스	
C19-162-965	30137805	Seven2Go 전극 클립	
C19-162-950	30122304	Seven2Go 손목 스트랩	
D08-162-930	51340291	InLab® 605 Polarographic DO Sensor	
D08-162-940	51344621	InLab® OptiOx DO Sensor	



S4-Standard Kit



S9-Standard Kit



30122300



51344621