

디지털 회전 농축기, RV10 Digital with Digital Oil Bath (HB Digital)

Cat. No.	Model	Glassware Type	EA Price ₩
E02-05-250	8041400	RV10 Digital V, 기본 수직형	10,770,000
E02-05-200	8041200	RV10 Digital Flex, without Glassware	6,850,000
E02-05-240	10004799	RV10 Digital V with 29/32 Glassware	7,600,000
E02-05-245	10004802	RV10 Digital V-C with 29/32 Glassware	8,150,000
E02-05-247	10006365	RV10 Digital with Dry Ice Condenser	8,670,000
E02-05-249	10006364	RV10 Digital with Dry Ice Condenser, Coated	9,820,000

- 탈취방식에 의한 장치의 편리한 탈부착
- 정방향, 역방향으로 Flask 회전 주기적 전환작동 가능
- Flask 회전속도 디지털 표시 및 타이머 기능 탑재
- Flask 회전속도 : 5 ~ 280 rpm (Digital Display)
- 앞면에 위치한 인체공학적 Operator를 통한 안전하고 편리한 조작
- Main Controller와 별도로 독립적으로 작동 가능하고, 다양한 안전회로로 이루어진 Heating Bath
- 불시에 전원이 Off 될 경우 Flask가 Heating Bath로부터 자동으로 상승하는 (140 mm) Safety Stop 기능
- Heating Bath 내 Flask의 위치를 자동인식하여 Flask의 깨짐 방지
- Bath의 손잡이가 내장된 Digital Oil Heating Bath
- Heating Bath의 온도가 Digital로 표시
- Heating Bath와 Drive Unit 사이의 적외선 데이터 전송방식
- Labworldsoft가 설치된 컴퓨터와 연결 가능한 RS232 Interface
- Labworldsoft 프로그램에 의한 자동 조작

진공 조절식 회전 농축기, RV10 Digital with Oil Bath (HB Digital) & Vacuum Controller & Vacuum Pump

Cat. No.	Model	Glassware Type	EA Price ₩
E02-05-204	10012017	RV10 Digital pro V with Pump, RC2 Lite Circulator	16,390,000
E02-05-207	10012018	RV10 Digital pro V with Pump, RC2 Green basic Circulator Double Complete	31,280,000

- VACSTAR Digital Pump
 - Chemical Resistant Membrane
 - 유량 22 L/min, 최대진공 2 mbar
 - Size : w150 x h375 x d370mm
 - Material in Contact with Medium : Al₂O₃, PTFE, FFPM, PPS, NBR

진공 조절식 회전 농축기, RV10 Auto Pro with Digital Oil Bath (HB Digital) & Vacuum Controller & Vacuum Pump

Cat. No.	Model	Glassware Type	EA Price ₩
E02-05-460	10000395	RV10 Auto Pro V with Pump	17,650,000
E02-05-410	10000397	RV10 Auto Pro FLEX with Pump, w/o Glassware	12,670,000
E02-05-470	10004815	RV 10 Auto Pro V with 29/32 Glassware	13,190,000
E02-05-475	10004812	RV10 Auto Pro V-C with 29/32 Glassware	14,240,000

- IKA VACSTAR Digital Pump는 제공된 데이터 케이블을 통해 RV 10 auto 메인 드라이브에서 직접 제어
- VACSTAR Digital Pump
 - Chemical Resistant Membrane
 - 유량 22 L/min, 최대진공 2 mbar
 - Size w150 x h375 x d370mm
 - Material in Contact with Medium : Al₂O₃, PTFE, FFPM, PPS, NBR



RV10 Digital



RV10 Auto Pro