

Mini Kuderna-Danish Evaporative Concentrator


미니 쿠데르나-데니스 농축기, LukeGL®

| Cat. No. | Model | Description | EA Price ₩ |
|-------------|---------|---------------|------------|
| E02-132-008 | L130450 | 40 ml 농축기 Set | 185,500 |

- For Concentration of Materials in Volatile Solvents
- Screw Type의 연결 방식으로 구성을 간소화 시킴 • LukeGL®

구성 Component

| Cat. No. | Model | Description | EA Price ₩ |
|-------------|---------|---------------------------------|------------|
| E02-132-607 | L130451 | Snyder Column 112 mm, 3 Section | 59,000 |
| E02-132-614 | L130452 | 40 ml Flask, 상 20-400, 하 13-425 | 37,100 |
| E02-132-621 | L130453 | 2 ml Receiver | 22,700 |
| E02-132-628 | L131823 | Connector, 20-400 to 20-400 | 64,900 |
| E02-132-635 | L130455 | Connector, 13-425 to 13-425 | 53,500 |



L130450

Kudema-Danish Evaporative Concentrator

쿠데르나-데니스 농축기, LukeGL®

| Cat. No. | Model | Description | Receiver | EA Price ₩ |
|-------------|---------|----------------|----------|------------|
| E02-132-016 | L130460 | 250 ml 농축기 Set | 10 ml | 108,300 |
| E02-132-024 | L130470 | 500 ml 농축기 Set | 15 ml | 112,500 |
| E02-132-032 | L130480 | 1 L 농축기 Set | 25 ml | 116,700 |

- Joint Size : 24/40 (Flask Top), 19/22 (Flask Bottom)
- Joint Clamp 24/40 및 19/22 용 각 2개씩 기본 포함 • LukeGL®

Spare Component

| Cat. No. | Model | Description | EA Price ₩ |
|-------------|---------|------------------------------------|------------|
| E02-132-701 | L130461 | Snyder Column 24/40 (공통) | 43,900 |
| E02-132-708 | L130462 | 250 ml Flask, 상 24/40, 하 19/22 | 37,100 |
| E02-132-715 | L130471 | 500 ml Flask, 상 24/40, 하 19/22 | 41,300 |
| E02-132-722 | L131881 | 1 L Flask, 상 24/40, 하 19/22 | 45,500 |
| E02-132-729 | L130463 | Receiver, 10 ml (Graduated, 19/22) | 20,600 |
| E02-132-736 | L130473 | Receiver, 15 ml (Graduated, 19/22) | 20,600 |
| E02-132-743 | L130483 | Receiver, 25 ml (Graduated, 19/22) | 20,600 |



쿠데르나-데니스 농축기

Receiver

Micro Soxhlet Extractor

미니 썩시렛 추출 장치, 25ml

| Cat. No. | Model | 용량 | Top Flange | Thimble | EA Price ₩ |
|-------------|---------|-------|------------|-----------|------------|
| E03-132-054 | 1374-01 | 25 ml | 19/38 | ∅10×50 mm | 554,700 |

- Thimble Filter (별도 구매) • Chemglass

구성 Component

| Cat. No. | Model | Description | EA Price ₩ |
|-------------|---------|---------------------------------|------------|
| E03-132-611 | 1374-02 | Extractor Body Only | 342,400 |
| E03-132-620 | 1374-03 | Condenser Only | 161,000 |
| F24-132-062 | 1506-86 | 25 ml, 14/20 Round Bottom Flask | 48,600 |

Optional Accessory

| Cat. No. | Model | 규격 | Size (mm) | Pore (um) | EA Price ₩ |
|-------------|---------|---------------|-----------|-----------|------------|
| E03-132-126 | 1374-04 | Glass Thimble | ∅10×h50 | 40 ~ 60 | 83,100 |
| E03-132-132 | 1374-05 | | | 170 ~ 220 | 86,700 |



Extractor, Liquid

액상 추출 장치

| Cat. No. | Model | 용량 | Stopcock | EA Price ₩ |
|-------------|-----------|-------------|----------|------------|
| E03-132-060 | 1360-01 ① | 1 ~ 1.2 L 용 | X | 823,000 |
| E03-132-066 | 1364-01 ② | | O | 1,013,700 |

- 물보다 무거운 액상용, EPA 법에 따름
- 상부 Outer Joint 45/50, 하부 Inner Joint 24/40 • Chemglass

Optional Accessory

| Cat. No. | Model | 규격 | EA Price ₩ |
|------------|---------|-----------------|------------|
| C18-10-462 | 52260G6 | 45/50 Condenser | 64,100 |

