

## Soxhlet Extractor for Solid 썩시렛 추출 장치, LukeGL®

Cat. No.	Model	용량	상부 Joint	Flask Joint	EA Price ₩
E03-132-006	L136801	100 ml	34/45	24/40	110,200
E03-132-012	L136802	250 ml	34/45		124,100
E03-132-018	L136803		45/50		162,700
E03-132-024	L136804	500 ml	34/45		134,900
E03-132-030	L136805		55/50		192,000
E03-132-036	L136806		1 L		71/60

- Borosilicate Glass 3.3
- Spare Round Bottom Flask는 452 Page 참조
- Thimble Filter (별도 구매)

### Extractor Only / 썩시렛 사이폰

Cat. No.	Model	Joint	용량	Thimble Model	EA Price ₩
E03-132-507	L136860	34/45	70 ml	2800288	48,900
E03-132-514	L136861	45/50	100 ml	2800338	71,800
E03-132-521	L136862	55/50	200 ml	2800432	100,100
E03-132-528	L136863	71/60	500 ml	10350273	219,300

### Allihn Condenser Only

Cat. No.	Model	Joint	길이	EA Price ₩
C18-10-450	L52260G4	34/45	170 mm	48,200
C18-10-462	L52260G6	45/50	200 mm	64,100
C18-10-750	L72260G4	55/50	250 mm	97,100
C18-10-756	L72260G6	71/60	300 mm	180,000

- Condenser의 길이는 냉각부만의 길이임

### Optional Thimble Filter, Cellulose

Cat. No.	Model	Size	적용 Joint	PK Price ₩
F15-07-205	2800288	Ø30×80 mm	34/45	25ea 302,000
F15-07-221	2800338	Ø35×80 mm	45/50	25ea 191,000
F15-07-237	2800432	Ø45×123 mm	55/50	25ea 234,000
F15-07-250	10350273	Ø52×145 mm	71/60	25ea 499,000

## Jumbo Soxhlet Extractor 대용량 썩시렛 추출 장치, Set

Cat. No.	Model	용량	Top Flange	Thimble Size	EA Price ₩
E03-132-038	1373-01	1 L	60 mm	Ø58×170 mm	2,064,100
E03-132-042	1373-02	3 L	100 mm	Ø81×230 mm	3,219,800
E03-132-048	1373-03	12 L	150 mm	Ø125×310 mm	5,034,000

- Borosilicate Glass 3.3
- O-ring형 Flange로 구성된 대용량 추출기
- Teflon O-ring 기본 포함
- Thimble Filter (별도 구매)
- Chemglass

Model	1373-01	1373-02	1373-03
Extractor Capacity	500 ml	1.5 L	5 L
Extractor id	68 mm	95 mm	140 mm
Extractor Flange	60	100	150
Lower joint	45/50	45/50	55/50
Condenser size	350 mm	400 mm	750 mm
Paper Thimble Size	Ø58×170 mm	Ø81×230 mm	Ø125×310 mm
Flask Capacity	1 L	3 L	12 L

### Optional Thimble Filter, Cellulose

Cat. No.	Model	Size	PK Price ₩
E03-132-050	1037-TF-01	Ø48×200 mm	25ea 581,000
E03-132-052	1037-TF-02	Ø75×250 mm	25ea 782,000
E03-132-056	1037-TF-03	Ø125×315 mm	1ea 1,103,800

