

Qualitative Filter Paper 정성여과지, Cellulose Fiber

- 정성분석, 생물학 수용액, 유기 용매 분석용
- 공기 모니터링, 식품분석
- GVS Lifescience



Grade	Pore Size	Thickness	Properties
1	10 ~ 13 um	180 um	medium
2	7 ~ 8 um	190 um	medium slow
3	5 ~ 7 um	320 um	medium-slow/thick
4	15 ~ 20 um	210 um	Very Fast
5	3 ~ 5 um	170 um	Very Slow

Grade 1 : Pore 10 ~ 13 um, 두께 180 um

Cat. No.	Model	Size	PK Price ₩
F12-594-004	FP01-042	Ø42.5mm	100ea 12,400
F12-594-008	FP01-047	Ø47 mm	100ea 12,400
F12-594-012	FP01-055	Ø55 mm	100ea 12,400
F12-594-016	FP01-070	Ø70 mm	100ea 13,500
F12-594-020	FP01-090	Ø90 mm	100ea 13,500
F12-594-024	FP01-110	Ø110 mm	100ea 17,500
F12-594-028	FP01-125	Ø125 mm	100ea 21,600
F12-594-032	FP01-150	Ø150 mm	100ea 30,900
F12-594-036	FP01-185	Ø185 mm	100ea 38,100

Grade 2 : Pore 7 ~ 8 um, 두께 190 um

Cat. No.	Model	Size	PK Price ₩
F12-594-058	FP02-042	Ø42.5mm	100ea 12,400
F12-594-060	FP02-047	Ø47 mm	100ea 12,400
F12-594-064	FP02-055	Ø55 mm	100ea 12,400
F12-594-068	FP02-070	Ø70 mm	100ea 13,500
F12-594-072	FP02-090	Ø90 mm	100ea 13,500
F12-594-076	FP02-110	Ø110 mm	100ea 17,500
F12-594-080	FP02-125	Ø125 mm	100ea 21,600
F12-594-084	FP02-150	Ø150 mm	100ea 30,900
F12-594-088	FP02-185	Ø185 mm	100ea 38,100

Grade 3 : Pore 5 ~ 7 um, 두께 320 um

Cat. No.	Model	Size	PK Price ₩
F12-594-090	FP03-042	Ø42.5 mm	100ea 19,600
F12-594-093	FP03-047	Ø47 mm	100ea 19,600
F12-594-096	FP03-055	Ø55 mm	100ea 19,600
F12-594-099	FP03-070	Ø70 mm	100ea 22,700
F12-594-102	FP03-090	Ø90 mm	100ea 24,700
F12-594-105	FP03-110	Ø110 mm	100ea 27,800
F12-594-108	FP03-125	Ø125 mm	100ea 34,000
F12-594-111	FP03-150	Ø150 mm	100ea 47,400
F12-594-114	FP03-185	Ø185 mm	100ea 73,100

Grade 4 : Pore 15 ~ 20 um, 두께 210 um

Cat. No.	Model	Size	PK Price ₩
F12-594-156	FP04-042	Ø42.5 mm	100ea 12,400
F12-594-160	FP04-047	Ø47 mm	100ea 12,400
F12-594-164	FP04-055	Ø55 mm	100ea 12,400
F12-594-168	FP04-070	Ø70 mm	100ea 13,500
F12-594-172	FP04-090	Ø90 mm	100ea 13,500
F12-594-176	FP04-110	Ø110 mm	100ea 17,500
F12-594-180	FP04-125	Ø125 mm	100ea 21,600
F12-594-184	FP04-150	Ø150 mm	100ea 30,900
F12-594-188	FP04-185	Ø185 mm	100ea 38,100

Grade 5 : Pore 3 ~ 5 um, 두께 170 um

Cat. No.	Model	Size	PK Price ₩
F12-594-200	FP05-042	Ø42.5 mm	100ea 12,400
F12-594-204	FP05-047	Ø47 mm	100ea 12,400
F12-594-208	FP05-055	Ø55 mm	100ea 12,400
F12-594-212	FP05-070	Ø70 mm	100ea 13,500
F12-594-216	FP05-090	Ø90 mm	100ea 13,500
F12-594-220	FP05-110	Ø110 mm	100ea 17,500
F12-594-224	FP05-125	Ø125 mm	100ea 21,600
F12-594-228	FP05-150	Ø150 mm	100ea 30,900
F12-594-232	FP05-185	Ø185 mm	100ea 38,100

Quantitative Filter Paper 정량여과지, Cellulose Fiber

- 시료 중에 함유된 여러 물질의 중량 측정 및 분석 여과용
- GVS Lifescience



Grade	Pore Size	Thickness	Properties
40	7 ~ 9 um	170 um	Medium slow
41	20 ~ 25 um	190 um	Fast
42	2 ~ 3 um	160 um	Very slow
44	2 ~ 4 um	160 um	Slow

Grade 40 : Pore 7 ~ 9 um, 두께 170 um

Cat. No.	Model	Size	PK Price ₩
F12-594-280	FP40-042	Ø42.5 mm	100ea 20,600
F12-594-284	FP40-047	Ø47 mm	100ea 20,600
F12-594-288	FP40-055	Ø55 mm	100ea 20,600
F12-594-292	FP40-070	Ø80 mm	100ea 23,700
F12-594-296	FP40-090	Ø90 mm	100ea 26,800
F12-594-300	FP40-110	Ø110 mm	100ea 28,800
F12-594-304	FP40-125	Ø125 mm	100ea 41,200
F12-594-308	FP40-150	Ø150 mm	100ea 57,700
F12-594-312	FP40-185	Ø185 mm	100ea 88,600

- 일반적 정량분석에 사용되는 Ashless Filter Paper

Grade 41 : Pore 20 ~ 25 um, 두께 190 um

Cat. No.	Model	Size	PK Price ₩
F12-594-360	FP41-042	Ø42.5 mm	100ea 20,600
F12-594-364	FP41-047	Ø47 mm	100ea 20,600
F12-594-368	FP41-055	Ø55 mm	100ea 20,600
F12-594-372	FP41-070	Ø80 mm	100ea 23,700
F12-594-376	FP41-090	Ø90 mm	100ea 26,800
F12-594-380	FP41-110	Ø110 mm	100ea 28,800
F12-594-384	FP41-125	Ø125 mm	100ea 41,200
F12-594-388	FP41-150	Ø150 mm	100ea 57,700
F12-594-392	FP41-185	Ø185 mm	100ea 88,600

- 가장 빠른 Ashless filter Paper로 입자가 큰 물질을 필터 시 효과적

다음 페이지에서 계속...