

Qualitative Filter Paper 정성여과지, Cellulose Fiber

- 분석 시료속에 존재하는 물성의 검출에 사용
- 실험실이나 생산현장에서 일반적인 여과용
- Ash Content : 0.1% 이하
- HD Micro

Grade	Description	Pore Size	Thickness
10	Medium Flow	6 ~ 10 μ m	170 μ m
20	Medium Flow	5 ~ 8 μ m	200 μ m
21	Medium Fast	8 ~ 12 μ m	200 μ m
22	Very Fast	12 ~ 15 μ m	210 μ m
100	Slow	3 ~ 5 μ m	160 μ m
200	Very Slow	2 ~ 3 μ m	160 μ m



Grade 10 : Pore 6 ~ 10 μ m, 두께 170 μ m

Cat. No.	Model	Size	PK Price ₩
F12-61-010	01006055	Ø55 mm	100ea 4,000
F12-61-014	01006070	Ø70 mm	100ea 4,300
F12-61-018	01006090	Ø90 mm	100ea 5,000
F12-61-022	01006110	Ø110 mm	100ea 5,800
F12-61-026	01006125	Ø125 mm	100ea 6,800
F12-61-030	01006150	Ø150 mm	100ea 8,400
F12-61-034	01006185	Ø185 mm	100ea 13,200
F12-61-038	01006200	Ø200 mm	100ea 18,800
F12-61-040	01006240	Ø240 mm	100ea 20,800
F12-61-042	01006250	Ø250 mm	100ea 26,600
F12-61-046	01006285	Ø285 mm	100ea 32,500
F12-61-050	01006300	Ø300 mm	100ea 36,600
F12-61-054	01006350	Ø350 mm	100ea 44,600

Grade 22 : Pore 12 ~ 15 μ m, 두께 210 μ m

Cat. No.	Model	Size	PK Price ₩
F12-61-154	02212055	Ø55 mm	100ea 5,700
F12-61-158	02212070	Ø70 mm	100ea 7,900
F12-61-162	02212090	Ø90 mm	100ea 9,300
F12-61-164	02212110	Ø110 mm	100ea 10,200
F12-61-168	02212125	Ø125 mm	100ea 13,200
F12-61-172	02212150	Ø150 mm	100ea 19,400
F12-61-176	02212185	Ø185 mm	100ea 26,900
F12-61-180	02212200	Ø200 mm	100ea 39,300
F12-61-184	02212240	Ø240 mm	100ea 43,300
F12-61-188	02212250	Ø250 mm	100ea 49,900
F12-61-192	02212285	Ø285 mm	100ea 57,500
F12-61-196	02212300	Ø300 mm	100ea 66,500
F12-61-200	02212350	Ø350 mm	100ea 87,300

Grade 20 : Pore 5 ~8 μ m, 두께 200 μ m

Cat. No.	Model	Size	PK Price ₩
F12-61-058	02005055	Ø55 mm	100ea 4,400
F12-61-062	02005070	Ø70 mm	100ea 4,700
F12-61-066	02005090	Ø90 mm	100ea 5,200
F12-61-070	02005110	Ø110 mm	100ea 6,100
F12-61-074	02005125	Ø125 mm	100ea 7,000
F12-61-078	02005150	Ø150 mm	100ea 9,100
F12-61-082	02005185	Ø185 mm	100ea 13,300
F12-61-086	02005200	Ø200 mm	100ea 21,700
F12-61-090	02005240	Ø240 mm	100ea 22,900
F12-61-094	02005250	Ø250 mm	100ea 29,300
F12-61-098	02005285	Ø285 mm	100ea 35,900
F12-61-102	02005300	Ø300 mm	100ea 38,300
F12-61-106	02005350	Ø350 mm	100ea 55,900

Grade 100 : Pore 3 ~ 5 μ m, 두께 160 μ m

Cat. No.	Model	Size	PK Price ₩
F12-61-204	10003055	Ø55 mm	100ea 5,700
F12-61-208	10003070	Ø70 mm	100ea 7,900
F12-61-212	10003090	Ø90 mm	100ea 9,300
F12-61-216	10003110	Ø110 mm	100ea 10,200
F12-61-220	10003125	Ø125 mm	100ea 13,200
F12-61-224	10003150	Ø150 mm	100ea 19,400
F12-61-228	10003185	Ø185 mm	100ea 26,900
F12-61-232	10003200	Ø200 mm	100ea 39,300
F12-61-236	10003240	Ø240 mm	100ea 43,300
F12-61-240	10003250	Ø250 mm	100ea 49,900
F12-61-244	10003285	Ø285 mm	100ea 57,500
F12-61-248	10003300	Ø300 mm	100ea 66,500
F12-61-252	10003350	Ø350 mm	100ea 87,300

Grade 21 : Pore 8 ~ 12 μ m, 두께 200 μ m

Cat. No.	Model	Size	PK Price ₩
F12-61-110	02110055	Ø55 mm	100ea 5,700
F12-61-114	02110070	Ø70 mm	100ea 7,900
F12-61-118	02110090	Ø90 mm	100ea 9,300
F12-61-122	02110110	Ø110 mm	100ea 10,200
F12-61-126	02110125	Ø125 mm	100ea 13,200
F12-61-130	02110150	Ø150 mm	100ea 19,400
F12-61-134	02110185	Ø185 mm	100ea 26,900
F12-61-138	02110200	Ø200 mm	100ea 39,300
F12-61-142	02110240	Ø240 mm	100ea 43,300
F12-61-146	02110250	Ø250 mm	100ea 49,900
F12-61-150	02110285	Ø285 mm	100ea 57,500

Grade 200 : Pore 2 ~ 3 μ m, 두께 160 μ m

Cat. No.	Model	Size	PK Price ₩
F12-61-256	20002055	Ø55 mm	100ea 5,700
F12-61-260	20002070	Ø70 mm	100ea 7,900
F12-61-264	20002090	Ø90 mm	100ea 9,300
F12-61-268	20002110	Ø110 mm	100ea 10,200
F12-61-272	20002125	Ø125 mm	100ea 13,200
F12-61-276	20002150	Ø150 mm	100ea 19,400
F12-61-280	20002185	Ø185 mm	100ea 26,900
F12-61-284	20002200	Ø200 mm	100ea 39,300
F12-61-288	20002240	Ø240 mm	100ea 43,300
F12-61-292	20002250	Ø250 mm	100ea 49,900
F12-61-296	20002285	Ø285 mm	100ea 57,500