

## Quantitative Filter Paper 정량여과지, Cellulose Fiber

- 시료중에 함유된 여러 물질의 중량 측정 및 분석 여과용
- Ash Content : 0.1% 이하
- HD Micro

Grade	Description	Pore Size	Thickness
50	Very Fast	12 ~ 15 $\mu\text{m}$	200 $\mu\text{m}$
51	Fast	8 ~ 12 $\mu\text{m}$	190 $\mu\text{m}$
52	Medium Slow	3 ~ 5 $\mu\text{m}$	160 $\mu\text{m}$
53	Very Slow	1 ~ 2 $\mu\text{m}$	200 $\mu\text{m}$
60	Slow	2 ~ 3 $\mu\text{m}$	150 $\mu\text{m}$
30	Medium Fast	5 ~ 8 $\mu\text{m}$	170 $\mu\text{m}$



### Grade 50 : Pore 12 ~ 15 $\mu\text{m}$ , 두께 200 $\mu\text{m}$

Cat. No.	Model	Size	PK Price ₩
F12-61-510	05012055	Ø55 mm	100ea 9,900
F12-61-514	05012070	Ø70 mm	100ea 11,600
F12-61-518	05012090	Ø90 mm	100ea 12,700
F12-61-522	05012110	Ø110 mm	100ea 18,000
F12-61-526	05012125	Ø125 mm	100ea 22,100
F12-61-530	05012150	Ø150 mm	100ea 27,800
F12-61-534	05012185	Ø185 mm	100ea 42,200
F12-61-538	05012200	Ø200 mm	100ea 74,700
F12-61-540	05012240	Ø240 mm	100ea 77,200
F12-61-542	05012250	Ø250 mm	100ea 97,600
F12-61-546	05012285	Ø285 mm	100ea 104,900
F12-61-550	05012300	Ø300 mm	100ea 124,100
F12-61-554	05012350	Ø350 mm	100ea 195,700

### Grade 53 : Pore 1 ~ 2 $\mu\text{m}$ , 두께 200 $\mu\text{m}$

Cat. No.	Model	Size	PK Price ₩
F12-61-654	05301055	Ø55 mm	100ea 9,900
F12-61-658	05301070	Ø70 mm	100ea 11,600
F12-61-662	05301090	Ø90 mm	100ea 13,300
F12-61-664	05301110	Ø110 mm	100ea 18,000
F12-61-668	05301125	Ø125 mm	100ea 22,100
F12-61-672	05301150	Ø150 mm	100ea 29,200
F12-61-676	05301185	Ø185 mm	100ea 44,500
F12-61-680	05301200	Ø200 mm	100ea 78,600
F12-61-684	05301240	Ø240 mm	100ea 81,100
F12-61-688	05301250	Ø250 mm	100ea 102,800
F12-61-692	05301285	Ø285 mm	100ea 110,400
F12-61-696	05301300	Ø300 mm	100ea 124,300
F12-61-700	05301350	Ø350 mm	100ea 195,700

### Grade 51 : Pore 8 ~ 12 $\mu\text{m}$ , 두께 190 $\mu\text{m}$

Cat. No.	Model	Size	PK Price ₩
F12-61-558	05108055	Ø55 mm	100ea 9,900
F12-61-562	05108070	Ø70 mm	100ea 11,600
F12-61-566	05108090	Ø90 mm	100ea 12,700
F12-61-570	05108110	Ø110 mm	100ea 18,000
F12-61-574	05108125	Ø125 mm	100ea 22,100
F12-61-578	05108150	Ø150 mm	100ea 27,800
F12-61-582	05108185	Ø185 mm	100ea 42,200
F12-61-586	05108200	Ø200 mm	100ea 74,700
F12-61-590	05108240	Ø240 mm	100ea 77,200
F12-61-594	05108250	Ø250 mm	100ea 97,600
F12-61-598	05108285	Ø285 mm	100ea 104,900
F12-61-602	05108300	Ø300 mm	100ea 124,100
F12-61-606	05108350	Ø350 mm	100ea 195,700

### Grade 60 : Pore 2 ~ 3 $\mu\text{m}$ , 두께 150 $\mu\text{m}$

Cat. No.	Model	Size	PK Price ₩
F12-61-704	06003055	Ø55 mm	100ea 9,900
F12-61-708	06003070	Ø70 mm	100ea 11,600
F12-61-712	06003090	Ø90 mm	100ea 12,700
F12-61-716	06003110	Ø110 mm	100ea 18,000
F12-61-720	06003125	Ø125 mm	100ea 22,100
F12-61-724	06003150	Ø150 mm	100ea 27,800
F12-61-728	06003185	Ø185 mm	100ea 42,200
F12-61-732	06003200	Ø200 mm	100ea 74,700
F12-61-736	06003240	Ø240 mm	100ea 77,200
F12-61-740	06003250	Ø250 mm	100ea 97,600
F12-61-744	06003285	Ø285 mm	100ea 104,900
F12-61-748	06003300	Ø300 mm	100ea 124,100
F12-61-752	06003350	Ø350 mm	100ea 195,700

### Grade 52 : Pore 3 ~ 5 $\mu\text{m}$ , 두께 160 $\mu\text{m}$

Cat. No.	Model	Size	PK Price ₩
F12-61-610	05204055	Ø55 mm	100ea 9,900
F12-61-614	05204070	Ø70 mm	100ea 11,600
F12-61-618	05204090	Ø90 mm	100ea 12,700
F12-61-622	05204110	Ø110 mm	100ea 18,000
F12-61-626	05204125	Ø125 mm	100ea 22,100
F12-61-630	05204150	Ø150 mm	100ea 27,800
F12-61-634	05204185	Ø185 mm	100ea 42,200
F12-61-638	05204200	Ø200 mm	100ea 74,700
F12-61-642	05204240	Ø240 mm	100ea 77,200
F12-61-646	05204250	Ø250 mm	100ea 97,600
F12-61-650	05204285	Ø285 mm	100ea 104,900

### Grade 30 : Pore 5 ~ 8 $\mu\text{m}$ , 두께 170 $\mu\text{m}$

Cat. No.	Model	Size	PK Price ₩
F12-61-756	03005055	Ø55 mm	100ea 9,900
F12-61-760	03005070	Ø70 mm	100ea 11,600
F12-61-764	03005090	Ø90 mm	100ea 12,700
F12-61-768	03005110	Ø110 mm	100ea 18,000
F12-61-772	03005125	Ø125 mm	100ea 22,100
F12-61-776	03005150	Ø150 mm	100ea 27,800
F12-61-780	03005185	Ø185 mm	100ea 42,200
F12-61-784	03005200	Ø200 mm	100ea 74,700
F12-61-788	03005240	Ø240 mm	100ea 77,200
F12-61-792	03005250	Ø250 mm	100ea 97,600
F12-61-796	03005285	Ø285 mm	100ea 104,900