

Qualitative Filter Paper 정성여과지, Cellulose Fiber

- 분석 시료속에 존재하는 물성의 검출에 사용
- 실험실이나 생산현장에서 일반적인 여과용
- Whatman

Grade	Description	Pore Size	Thickness
1	Medium Flow	11 um	180 um
2	Medium Flow	8 um	190 um
3	Medium Flow, Thick	6 um	390 um
4	Very Fast	20 ~ 25 um	205 um
5	Slow	2.5 um	200 um
6	Medium to Slow	3 um	180 um

- Particle Retention Rating at 98% Efficiency



Grade 1 : Pore 11 um, 두께 180 um

Cat. No.	Model	Size	PK Price ₩
F12-07-004	1001325	Ø25 mm	100ea 19,600
F12-07-008	1001047	Ø47 mm	100ea 18,900
F12-07-012	1001055	Ø55 mm	100ea 17,100
F12-07-016	1001070	Ø70 mm	100ea 17,100
F12-07-020	1001090	Ø90 mm	100ea 12,800
F12-07-024	1001110	Ø110 mm	100ea 13,000
F12-07-028	1001125	Ø125 mm	100ea 23,100
F12-07-032	1001150	Ø150 mm	100ea 22,900
F12-07-036	1001185	Ø185 mm	100ea 43,500
F12-07-040	1001240	Ø240 mm	100ea 82,100
F12-07-044	1001270	Ø270 mm	100ea 98,600
F12-07-048	1001320	Ø320 mm	100ea 125,000
F12-07-052	1001929	w600×d600 mm	100ea 376,000

Grade 4 : Pore 20 ~ 25 um, 두께 205 um

Cat. No.	Model	Size	PK Price ₩
F12-07-156	1004047	Ø47 mm	100ea 26,200
F12-07-160	1004055	Ø55 mm	100ea 17,100
F12-07-164	1004070	Ø70 mm	100ea 17,100
F12-07-168	1004090	Ø90 mm	100ea 17,100
F12-07-172	1004110	Ø110 mm	100ea 17,800
F12-07-176	1004125	Ø125 mm	100ea 29,800
F12-07-180	1004150	Ø150 mm	100ea 35,900
F12-07-184	1004185	Ø185 mm	100ea 52,600
F12-07-188	1004240	Ø240 mm	100ea 102,000
F12-07-192	1004270	Ø270 mm	100ea 248,000
F12-07-196	1004320	Ø320 mm	100ea 357,000

Grade 2 : Pore 8 um, 두께 190 um

Cat. No.	Model	Size	PK Price ₩
F12-07-058	1002042	Ø42.5 mm	100ea 17,100
F12-07-060	1002047	Ø47 mm	100ea 22,300
F12-07-064	1002055	Ø55 mm	100ea 17,800
F12-07-068	1002070	Ø70 mm	100ea 18,900
F12-07-072	1002090	Ø90 mm	100ea 14,900
F12-07-076	1002110	Ø110 mm	100ea 14,700
F12-07-080	1002125	Ø125 mm	100ea 26,600
F12-07-084	1002150	Ø150 mm	100ea 24,200
F12-07-088	1002185	Ø185 mm	100ea 34,300
F12-07-092	1002240	Ø240 mm	100ea 93,300
F12-07-096	1002270	Ø270 mm	100ea 212,000
F12-07-100	1002320	Ø320 mm	100ea 216,000
F12-07-104	1002929	w600×d600 mm	100ea 494,000

Grade 5 : Pore 2.5 um, 두께 200 um

Cat. No.	Model	Size	PK Price ₩
F12-07-200	1005325	Ø25 mm	100ea 35,700
F12-07-204	1005047	Ø47 mm	100ea 23,200
F12-07-208	1005055	Ø55 mm	100ea 21,400
F12-07-212	1005070	Ø70 mm	100ea 23,200
F12-07-216	1005090	Ø90 mm	100ea 23,400
F12-07-220	1005110	Ø110 mm	100ea 26,700
F12-07-224	1005125	Ø125 mm	100ea 40,600
F12-07-228	1005150	Ø150 mm	100ea 50,100
F12-07-232	1005185	Ø185 mm	100ea 84,100
F12-07-236	1005240	Ø240 mm	100ea 257,000
F12-07-244	1005320	Ø320 mm	100ea 441,000

Grade 3 : Pore 6 um, 두께 390 um

Cat. No.	Model	Size	PK Price ₩
F12-07-112	1003055	Ø55 mm	100ea 29,000
F12-07-116	1003070	Ø70 mm	100ea 32,800
F12-07-120	1003090	Ø90 mm	100ea 35,300
F12-07-124	1003110	Ø110 mm	100ea 40,000
F12-07-128	1003125	Ø125 mm	100ea 47,100
F12-07-132	1003150	Ø150 mm	100ea 63,600
F12-07-136	1003185	Ø185 mm	100ea 87,300
F12-07-140	1003240	Ø240 mm	100ea 136,000
F12-07-148	1003320	Ø320 mm	100ea 515,000

Grade 6 : Pore 3 um, 두께 180 um

Cat. No.	Model	Size	PK Price ₩
F12-07-248	1006042	Ø42.5 mm	100ea 19,900
F12-07-252	1006070	Ø70 mm	100ea 23,100
F12-07-256	1006090	Ø90 mm	100ea 24,800
F12-07-260	1006110	Ø110 mm	100ea 23,200
F12-07-264	1006125	Ø125 mm	100ea 37,300
F12-07-268	1006150	Ø150 mm	100ea 48,400
F12-07-272	1006185	Ø185 mm	100ea 86,900
F12-07-276	1006240	Ø240 mm	100ea 143,000