

Quantitative Filter Paper 정량여과지, Cellulose Fiber

- 시료중에 함유된 여러 물질의 중량 측정 및 분석 여과용
- Whatman

Grade	Description	Pore Size	Thickness
40	Medium	8 um	210 um
41	Fast	20 ~ 25 um	220 um
42	Slow	2.5 um	200 um
43	Medium to Fast	16 um	220 um
44	Slow to Medium	3 um	180 um

- Ash Content : 0.007% (at 900°C in Air)
- Particle Retention Rating at 98% Efficiency

Grade 40 : Pore 8 um, 두께 210 um

Cat. No.	Model	Size	PK Price ₩
F12-07-280	1440047	Ø47 mm	100ea 47,800
F12-07-284	1440055	Ø55 mm	100ea 32,300
F12-07-288	1440070	Ø70 mm	100ea 33,900
F12-07-292	1440090	Ø90 mm	100ea 42,700
F12-07-296	1440110	Ø110 mm	100ea 40,600
F12-07-300	1440125	Ø125 mm	100ea 61,600
F12-07-304	1440150	Ø150 mm	100ea 82,400
F12-07-308	1440185	Ø185 mm	100ea 150,000
F12-07-312	1440240	Ø240 mm	100ea 258,000
F12-07-316	1440320	Ø320 mm	100ea 523,000

Grade 41 : Pore 20 ~ 25 um, 두께 220 um

Cat. No.	Model	Size	PK Price ₩
F12-07-320	1441047	Ø47 mm	100ea 31,400
F12-07-324	1441055	Ø55 mm	100ea 33,900
F12-07-328	1441070	Ø70 mm	100ea 33,900
F12-07-332	1441090	Ø90 mm	100ea 44,100
F12-07-336	1441110	Ø110 mm	100ea 44,500
F12-07-340	1441125	Ø125 mm	100ea 70,100
F12-07-344	1441150	Ø150 mm	100ea 86,900
F12-07-348	1441185	Ø185 mm	100ea 131,000
F12-07-352	1441240	Ø240 mm	100ea 239,000
F12-07-356	1441320	Ø320 mm	100ea 435,000

Grade 42 : Pore 2.5 um, 두께 200 um

Cat. No.	Model	Size	PK Price ₩
F12-07-360	1442047	Ø47 mm	100ea 44,400
F12-07-364	1442055	Ø55 mm	100ea 30,700
F12-07-368	1442070	Ø70 mm	100ea 35,100
F12-07-372	1442090	Ø90 mm	100ea 44,400
F12-07-376	1442110	Ø110 mm	100ea 41,500
F12-07-380	1442125	Ø125 mm	100ea 63,200
F12-07-384	1442150	Ø150 mm	100ea 80,600
F12-07-388	1442185	Ø185 mm	100ea 131,000
F12-07-392	1442240	Ø240 mm	100ea 245,000
F12-07-396	1442320	Ø320 mm	100ea 503,000

Grade 43 : Pore 16 um, 두께 220 um

Cat. No.	Model	Size	PK Price ₩
F12-07-400	1443090	Ø90 mm	100ea 103,000
F12-07-404	1443110	Ø110 mm	100ea 81,800
F12-07-408	1443125	Ø125 mm	100ea 145,000
F12-07-412	1443150	Ø150 mm	100ea 183,000
F12-07-416	1443185	Ø185 mm	100ea 270,000



Grade 44 : Pore 3 um, 두께 180 um

Cat. No.	Model	Size	PK Price ₩
F12-07-420	1444070	Ø70 mm	100ea 73,900
F12-07-424	1444090	Ø90 mm	100ea 56,300
F12-07-428	1444110	Ø110 mm	100ea 57,800
F12-07-432	1444125	Ø125 mm	100ea 145,000
F12-07-436	1444150	Ø150 mm	100ea 124,000
F12-07-440	1444185	Ø185 mm	100ea 270,000
F12-07-444	1444240	Ø240 mm	100ea 428,000

Hardened Ashless Filter Paper

Grade	Description	Pore Size	Thickness
540	Medium	8 um	160 um
541	Fast	20 ~ 25 um	155 um
542	Slow	2.7 um	150 um

- Ash Content : ≤0.006%
- Particle Retention Rating at 98% Efficiency
- High Wet Strength와 Chemical Resistance

Grade 540 : Pore 8 um, 두께 160 um

Cat. No.	Model	Size	PK Price ₩
F12-07-448	1540055	Ø55 mm	100ea 65,000
F12-07-456	1540090	Ø90 mm	100ea 68,900
F12-07-460	1540110	Ø110 mm	100ea 91,100
F12-07-464	1540125	Ø125 mm	100ea 100,000
F12-07-468	1540150	Ø150 mm	100ea 211,000
F12-07-472	1540185	Ø185 mm	100ea 312,000

Grade 541 : Pore 22 um, 두께 155 um

Cat. No.	Model	Size	PK Price ₩
F12-07-476	1541055	Ø55 mm	100ea 43,300
F12-07-480	1541070	Ø70 mm	100ea 46,900
F12-07-488	1541090	Ø90 mm	100ea 107,000
F12-07-492	1541110	Ø110 mm	100ea 84,400
F12-07-496	1541125	Ø125 mm	100ea 99,400
F12-07-500	1541150	Ø150 mm	100ea 141,000
F12-07-504	1541185	Ø185 mm	100ea 312,000

Grade 542 : Pore 2.7 um, 두께 150 um

Cat. No.	Model	Size	PK Price ₩
F12-07-508	1542055	Ø55 mm	100ea 68,100
F12-07-512	1542070	Ø70 mm	100ea 81,200
F12-07-516	1542090	Ø90 mm	100ea 107,000
F12-07-520	1542110	Ø110 mm	100ea 87,800
F12-07-524	1542125	Ø125 mm	100ea 100,000
F12-07-528	1542150	Ø150 mm	100ea 211,000
F12-07-532	1542185	Ø185 mm	100ea 312,000