

Glass Microfiber Filter

유리 섬유 여지

- 재질 : 100% Borosilicate Glass
- Binder - free Glass Microfiber
- Max. Temperature : 500°C
- Whatman

Grade	Description	Pore Size	Thickness
GF/A	Fast, High Loading	1.6 um	260 um
GF/B	Medium to Fast, Very High Loading	1.0 um	675 um
GF/C	Medium to Fast, High Loading	1.2 um	260 um
GF/D	Fast, Very High Loading	2.7 um	675 um
GF/F	Medium, High Loading	0.7 um	420 um

- Particle Retention Rating at 98% efficiency



Grade GF/A : Pore 1.6 um, 두께 260 um

Cat. No.	Model	Size	PK Price ₩
F12-07-536	1820025	Ø25 mm	100ea 54,000
F12-07-540	1820047	Ø47 mm	100ea 66,300
F12-07-544	1820055	Ø55 mm	100ea 93,800
F12-07-548	1820070	Ø70 mm	100ea 110,000
F12-07-552	1820090	Ø90 mm	100ea 130,000
F12-07-554	1820110	Ø110 mm	100ea 148,000
F12-07-556	1820125	Ø125 mm	100ea 221,000
F12-07-560	1820150	Ø150 mm	100ea 280,000

Grade GF/D : Pore 2.7 um, 두께 675 um

Cat. No.	Model	Size	PK Price ₩
F12-07-620	1823025	Ø25 mm	100ea 60,100
F12-07-624	1823047	Ø47 mm	100ea 115,000
F12-07-628	1823055	Ø55 mm	100ea 180,000
F12-07-632	1823070	Ø70 mm	100ea 190,000
F12-07-636	1823090	Ø90 mm	25ea 103,000
F12-07-638	1823110	Ø110 mm	25ea 110,000
F12-07-640	1823125	Ø125 mm	25ea 180,000
F12-07-644	1823150	Ø150 mm	25ea 310,000

Grade GF/B : Pore 1.0 um, 두께 675 um

Cat. No.	Model	Size	PK Price ₩
F12-07-564	1821025	Ø25 mm	100ea 168,000
F12-07-568	1821047	Ø47 mm	100ea 103,000
F12-07-572	1821055	Ø55 mm	100ea 170,000
F12-07-576	1821070	Ø70 mm	100ea 200,000
F12-07-580	1821090	Ø90 mm	25ea 99,500
F12-07-582	1821110	Ø110 mm	25ea 115,000
F12-07-584	1821125	Ø125 mm	25ea 195,000
F12-07-588	1821150	Ø150 mm	25ea 184,000

Grade GF/F : Pore 0.7 um, 두께 420 um

Cat. No.	Model	Size	PK Price ₩
F12-07-648	1825025	Ø25 mm	100ea 102,000
F12-07-652	1825047	Ø47 mm	100ea 141,000
F12-07-656	1825055	Ø55 mm	100ea 482,000
F12-07-660	1825070	Ø70 mm	100ea 334,000
F12-07-664	1825090	Ø90 mm	25ea 172,000
F12-07-666	1825110	Ø110 mm	25ea 181,000
F12-07-668	1825125	Ø125 mm	25ea 348,000
F12-07-672	1825150	Ø150 mm	25ea 339,000

Grade GF/C : Pore 1.2 um, 두께 260 um

Cat. No.	Model	Size	PK Price ₩
F12-07-592	1822025	Ø25 mm	100ea 47,900
F12-07-596	1822047	Ø47 mm	100ea 39,800
F12-07-600	1822055	Ø55 mm	100ea 108,000
F12-07-604	1822070	Ø70 mm	100ea 144,000
F12-07-608	1822090	Ø90 mm	100ea 145,000
F12-07-610	1822110	Ø110 mm	100ea 119,000
F12-07-612	1822125	Ø125 mm	100ea 259,000
F12-07-616	1822150	Ø150 mm	100ea 362,000