

PP Capsule Filter, Capsflow

PP 캡슐 필터

Cat. No.	Model	Pore	Size	Fitting In/Out	EA Price ₩
F16-594-003	CIKPP10GLL5H	1.0 um	5"	1.5" Tri-clamp	85,000
F16-594-006	CIKPP30GLL5H	3.0 um			85,000
F16-594-009	CIKPP50GLL5H	5.0 um			85,000
F16-594-033	CSKPP02G051N	0.2 um	1.5" (500cm ² filtration area)	1" Tri-clamp	48,000
F16-594-039	CSKPP10G051N	1 um			51,000
F16-594-045	CSKPP30G051N	3 um	95% filtering	1"	51,000
F16-594-051	CSKPP50G051N	5 um			51,000
F16-594-036	CSKPP02P051N	0.2 um			1.5" (500cm ² filtration area)
F16-594-042	CSKPP10P051N	1 um	59,000		
F16-594-048	CSKPP30P051N	3 um	99% filtering	1"	59,000
F16-594-054	CSKPP50P051N	5 um			59,000

- 화학적 호환성이 광범위
- 내구성이 좋고 필터링된 잔여물 보유 용량이 큼
- 식초, 수용액, 맥주, 와인 및 증류주, 식음료, 식용유
- 수(水)공정, 화학물질 대량생산, 제약 중간제
- GVS Lifescience



PES Capsule Filter, Capsflow

PES 캡슐 필터

Cat. No.	Model	Pore	Size	Fitting In/Out	EA Price ₩
F16-594-013	CIKPS02BLL5H	0.2 um	5"	1.5" Tri-clamp	172,000
F16-594-016	CIKPS04BLL5H	0.45 um			172,000
F16-594-063	CSKPS02B051N	0.2 um	1.5" (Low Bio-burden)	1" Tri-clamp	96,000
F16-594-069	CSKPS04B051N	0.45 um			96,000
F16-594-066	CSKPS02G051N	0.2 um	1.5" (500cm ² filtration area)	1"	82,000
F16-594-072	CSKPS04G051N	0.45 um			82,000
F16-594-075	CSKPS40S051N	0.45 um + 0.2 um			135,000
F16-594-078	CSKPS80S051N	0.8 um + 0.45 um			135,000

- 매체에 흡착되는 물질이 포함된 여과에 적합
- 백신 및 생물학적 물질과 같은 단백질 용액 여과에 적합
- 의약품(세포배양, 대용량 주사, 진단, 혈액-혈청 분출), 식음료 정제 등
- GVS Lifescience



PTFE Capsule Filter, Capsflow

PTFE 캡슐 필터

Cat. No.	Model	Pore	Size	Fitting In/Out	EA Price ₩
F16-594-023	CIKPT02BLL5H	0.2 um	5"	1.5" Tri-clamp	165,000
F16-594-083	CSKPT02B051N	0.22 um	1.5" (500cm ² filtration area)	1" Tri-clamp	80,000
F16-594-086	CSKPT02G051N	0.22 um			75,000
F16-594-089	CSKPT04B051N	0.45 um			80,000
F16-594-092	CSKPT04G051N	0.45 um			75,000

- Hydrophobic, 내구성이 좋음
- 멸균 공기 공급, 의약품, 용매 등으로 사용
- FDA 규정 준수
- GVS Lifescience



Steaming in Place Sterilization Capsule Filter

SIP 멸균 캡슐 필터, Capsflow

Cat. No.	Model	Material	Pore	Fitting In/Out	EA Price ₩
F16-594-103	CXKNY01XLL5H	Nylon	0.1 um	1.5" Tri-clamp	418,000
F16-594-106	CXKNY02XLL5H		0.2 um		418,000
F16-594-113	CXKPS01XLL5H	PES	0.1 um		455,000
F16-594-116	CXKPS02XLL5H		0.2 um		455,000
F16-594-123	CXKPT00XLL5H	PTFE	0.05 um		484,000
F16-594-126	CXKPT01XLL5H		0.1 um		440,000
F16-594-129	CXKPT02XLL5H		0.2 um	440,000	

- ISO 9001 기준에 따라 제조
- Sanitary Flange 와 쉽게 연결 가능
- 액체 및 기체의 멸균 여과
- FDA 규정 준수
- Size : 5 inch
- GVS Lifescience

